


# FCI Nunchaku®

## Self-Retaining Nasolacrimal Duct Intubation Without Nasal Recovery

### INDICATIONS

- Canalicular Pathologies
- Dacryocystorhinostomy (DCR)
- Congenital Lacrimal Duct Obstruction



**FCI**®

a ZEISS company



[www.FCI-Ophthalmics.com](http://www.FCI-Ophthalmics.com)

## PRESENTATION

Nunchaku is a pushed silicone self-retaining bicanalicular intubation stent that acts like a conformer, allowing tears to be drained by capillarity, no retrieval from the nose is needed. The metallic guide introducer is located inside the lumen and not as an extension of the stent as in conventional intubation sets. The stability is guaranteed by the design of the silicone tubes, no need for knots or sutures in the nasal fossa.


<b>S1.1361</b>	90 mm (for children)	Box of 1, sterile
<b>S1.1371</b>	105 mm (for adults)	Box of 1, sterile

## CHARACTERISTICS

### NUNCHAKU TUBES


#### 10 MM MARK:

Distance between the punctum and the beginning of the lacrimal sac.




#### 15 MM MARK:

Distance between the punctum and the end of the lacrimal sac.



### METALLIC GUIDES:

The metallic guides give rigidity to the Nunchaku tubes and facilitate insertion in the lacrimal system.




## INITIAL PROBING

### DIAGNOSIS

Complex stenosis (contraindication) is distinguished from scarred nasolacrimal stenosis by tactile probing.

### DETECTING FALSE PASSAGES

A second, wider lacrimal probe with a blunt tip is introduced through the naris and gently steered through the inferior nasal meatus until metal-to-metal contact is achieved.



Probing

Metal-to-Metal Contact


### SELECTION OF STENT LENGTH

The selection of the stent length depends on the surgeon's preference. It is recommended to use a 90 mm stent for children and a 105 mm stent for adults in cases of classic intubation.

## IMPLANTING THE STENT


- 1 Stent is inserted into the canaliculus through the dilated punctum.
- 2 At bony contact, stent is rotated 90° and vertical catheterization is continued.
- 3 Once the nasal fossa floor is reached, the metallic guide is gently removed while the silicone tube is maintained in place.

The same procedure is repeated for the second canaliculus. A self-retaining bicanalicular intubation is achieved. No knots or sutures are needed at the end of the procedure.



## SPECIAL CASES

In certain special cases (as indicated to the right), remove the metallic guide a few millimeters. Then, push the silicone stent together with the metallic guide using forceps. This way, the silicone stent will not create a false passage.




Partial Stenosis of the Canaliculus

Major Curvature of the Nasolacrimal Duct

## INSERTION ERRORS

If the stent takes a wrong direction during the insertion, it risks making a false passage. To avoid this, the implantation procedure must be carefully followed.



False Passage in the Canaliculus

False Passage in the Nasolacrimal Duct



**FCI OPHTHALMICS**  
30 Corporate Park Drive, Suite 310/320  
Pembroke, MA 02359

**EMAIL**  
[info@fci-ophthalmics.com](mailto:info@fci-ophthalmics.com)

**TELEPHONE**  
800.932.4202 • 781.826.9060

**FAX**  
781.826.9062

**STAY CONNECTED!**



Social Media



Learn More!

[www.FCI-Ophthalmics.com](http://www.FCI-Ophthalmics.com)