

FCI Ophthalmics

PRODUCT CATALOG



- **Dry Eye**
- **Lacrimal Stents & Tubes**
- **Eyelid Repair**
- **Orbital**
- **Cataract**
- **Vitreoretinal**

FCI-Ophthalmics.com
800.932.4202



THE FCI STORY

Since its inception in 1984, FCI has become a world leader in products that solve problems in the areas of Dry Eye, Oculoplastics, Cataract and Retinal procedures.

Through close collaboration with leading ophthalmic surgeons worldwide, FCI is continually developing and testing products that contribute to the quality of patients' lives and products that respond to physicians' needs with advanced materials and innovative research.

FCI is excited about the future and is dedicated to the advancement of quality in eye care. FCI conforms to the world's standards for quality in design, manufacturing and service:

All FCI Ophthalmics products approved by FDA for commercial distribution.

FCI-Ophthalmics.com

800.932.4202

TABLE OF CONTENTS

Dry Eye

- 2 PERMANENT PUNCTUM PLUGS
- 4 TEMPORARY PUNCTUM PLUGS
- 5 TEAR TEST
- 5 DRY EYE ACCESSORIES

Lacrimal Stents & Tubes

- 6 BALLOON CATHETER
- 7 BICANALICULAR INTUBATION
- 9 MONOCANALICULAR INTUBATION
- 10 JONES TUBES
- 11 LACRIMAL STENTS & TUBES ACCESSORIES

Eyelid Repair

- 12 PTOSIS REPAIR
- 13 IMPLANTS FOR LAGOPHTHALMOS
- 14 EXTERNAL WEIGHTS FOR LAGOPHTHALMOS
- 15 EYELID REPAIR ACCESSORIES

Orbital

- 16 BIOCERAMIC IMPLANTS
- 16 SILICONE IMPLANTS
- 17 SELF-INFLATING TISSUE EXPANDERS
- 18 ORBITAL ACCESSORIES

Cataract

- 21 CAPSULAR TENSION RINGS
- 23 CATARACT ACCESSORIES
- 25 KITARO KITS

Vitreoretinal

- 26 TROCAR CANNULA SYSTEM
- 27 PROBES
- 30 BACKFLUSH INSTRUMENTS
- 31 VITREORETINAL INSTRUMENTS
- 33 SCLERAL BUCKLING COMPONENTS
- 34 VITREORETINAL ACCESSORIES
- 35 INTRAVITREAL INJECTION DEVICE
- 36 INDEX

One Plug One Size

SNUGPLUG™

That's *all* you need!

One-size-fits-all punctum plug **eliminates sizing** and **simplifies stocking**

Pre-stretched shape (on inserter) **avoids dilation** and **facilitates insertion**

Expanded shape (once inserted) assures snug fit and **virtually eliminates pop-out**

Soft collar prevents migration and **provides patient comfort**



Visit FCI-Ophthalmics.com to watch the SnugPlug informational video.



Permanent Punctum Plugs

SNUGPLUG™

Plugs are packaged preloaded in a stretched position and return to their natural shape when released from inserter.

- For punctal occlusion
- Medical grade silicone
- Permanent yet removable
- Individually packaged

2 preloaded plugs per box, sterile

S2.4002u One size fits all

FCI
EXCLUSIVE



Natural Shape of SnugPlug in Punctum



Pre-stretched SnugPlug

Permanent Punctum Plugs (cont.)

READY-SET PLUGS

Designed by J.A. Bernard, M.D.

Ultra-thin collarette conforms to natural anatomy of the eyelid for maximum patient comfort.

Preloaded on disposable inserter/dilator instrument for ease of insertion. Easily removed with forceps.

- For punctal occlusion
- Medical grade silicone
- 2 mm long
- Individually packaged

2 preloaded plugs per box, sterile

Slim Plugs

The Slim Mini and Slim Petite plugs offer a slightly smaller bulb for easier insertion into small puncta.

Exclusive Slanted Collarette Plugs

The Micro to Large plugs feature the unique slanted collarette that hugs the natural curve of the eyelid.


The slightly larger bulb resists migration and is designed to prevent rub-out.

SLIM PLUGS

S2.3632u	0.4 mm	Slim Mini
S2.3642u	0.5 mm	Slim Petite

SLANTED COLLARETTE PLUGS

S2.3652u	0.6 mm	Micro
S2.3672u	0.7 mm	Mini
S2.3422u	0.8 mm	Small
S2.3122u	0.9 mm	Medium
S2.3522u	1.0 mm	Large



NON-STERILE PLUGS

Ready-Set plugs are pre-mounted on a wheel for quick and easy loading onto a plug inserter.

- For punctal occlusion
- Medical grade silicone

24 plugs per box, 2 inserters, 1 wheel, non-sterile

Variety Pack: 6 plugs of each size (0.4, 0.5, 0.6, 0.7 mm)

SLIM PLUGS


S2.3324u	0.4 mm	Slim Mini
S2.3424u	0.5 mm	Slim Petite

SLANTED COLLARETTE PLUGS

S2.3524u	0.6 mm	Micro
S2.3724u	0.7 mm	Mini

VARIETY PACK

S2.3624u	Assorted	6 plugs each
-----------------	----------	--------------



TEARSAVER PLUGS

The tapered shaft design helps to anchor the plug securely in the punctum for more comfort and less risk of rub-out.

Preloaded on disposable inserter/dilator instrument for ease of insertion.

- For punctal occlusion
- Medical grade silicone
- Permanent yet removable
- Individually packaged

2 preloaded plugs per box, sterile

TAPERED PLUGS

814004	0.4 mm
814005	0.5 mm
814006	0.6 mm
814007	0.7 mm
814008	0.8 mm



Permanent Punctum Plugs (cont.)

PVP PERFORATED PLUGS

Designed by J.A. Bernard, M.D.

These plugs are coated with a thin layer of polyvinylpyrrolidone (PVP) that makes the surface slick to prevent debris from collecting on its surface and enhance drainage.

- For partial occlusion or punctal stenosis
- Medical grade silicone with PVP coating
- Individually packaged

2 preloaded plugs per box, sterile

MINI PLUG

S1.3522u 0.7 mm 0.4 Diameter Channel

MEDIUM PLUG

S1.3512u 0.9 mm 0.6 Diameter Channel



Perforation through entire plug allows tear drainage through blocked punctum.

Temporary Punctum Plugs

COLLAGEN PLUGS

- To diagnose dry eye
- To assess the effectiveness of permanent punctal occlusion
- Collagen-dissolvable in 5–7 days
- 2 mm long

6 plugs per packet/12 packets per box (72 total), sterile

Variety Pack: 6 plugs per packet/4 packets of each size per box (72 total), sterile

CP-0002 0.2 mm

CP-0003 0.3 mm

CP-0004 0.4 mm

CP-0234 Variety Pack



PROLONG™ ABSORBABLE LONG TERM PLUGS

- For short term dry eye conditions
- For post surgical dry eye symptoms
- Absorbable copolymer – lasts up to 3 months
- 2 mm long

2 plugs per packet, 10 packets per box (20 total), sterile

Variety Pack: 3 packets 0.3 mm, 5 packets 0.4 mm, 2 packets 0.5 mm, sterile

TP-3000 0.3 mm

TP-4000 0.4 mm

TP-5000 0.5 mm

TP-3450 Variety Pack



Tear Test

ZONE-QUICK

A simple and accurate system that measures tear volume in just 15 seconds per eye. Each box has an easy to read diagnostic indicator ruler.

- For tear volume measurement
- Phenol red thread tear test
- Results ready in 15 seconds

50 pairs per box, sterile

ZQ-1010 100 Threads



Dry Eye Accessories

PUNCTAL GAUGES FOR READY-SET PLUGS

Each instrument has one gauge size on one end, and a different gauge size on the opposite end.

- For proper plug sizing
- Stainless steel
- Includes autoclavable sterilizing case

Box of 1, non-sterile

PUNCTAL GAUGE SET

S2.3060u 6 Gauge Sizes (3 Instruments)




REUSABLE DOUBLE DILATOR

- Stainless steel

Box of 1, non-sterile

A8.4050 Fine & Medium Tapers



Lacrimal Duct Balloon Catheter

OPHTACATH®

Fast and Efficient Procedure

OphtaCath® achieves true dilation of the lacrimal duct and effectively treats the symptoms of epiphora.

- Eliminates placement of tubes
- Less invasive, reduced trauma



Visit FCI-Ophthalmics.com to watch the OphtaCath informational video.

Inflation Device

Balloon Inflated

Balloon Catheter

OPHTACATH®

OphtaCath® achieves true dilation of the lacrimal duct and effectively treats the symptoms of epiphora.

- For nasolacrimal duct obstructions
- Fast effective procedure
- Eliminates placement of tubes
- Less invasive, reduced trauma
- Bilateral includes stopcock for simultaneous inflation

Box of 1, sterile



2 MM and 3 MM UNILATERAL KIT

S1.4121 1x2 mm Balloon Catheter & Disposable Inflation Device

S1.4131 1x3 mm Balloon Catheter & Disposable Inflation Device

2 MM and 3 MM BILATERAL KIT

S1.4122 2x2 mm Balloon Catheters & Disposable Inflation Device

S1.4132 2x3 mm Balloon Catheters & Disposable Inflation Device

Bicanalicular Intubation

NUNCHAKU®

- A pushed self-retaining bicanalicular stent
- For canalicular pathologies
- For dacryocystorhinostomy (DCR)
- For nasal lacrimal duct obstruction
- Stability through silicone tube design; no knot or suture needed

Box of 1, sterile

S1.1361 90 mm (for Children)

S1.1371 105 mm (for Adults)

FCI
EXCLUSIVE



SELF-RETAINING SET II (SRS II)

Designed by P. Bigé, M.D.

- For horizontal lacrimal duct obstruction
- Punctal stenosis
- Silicone
- Set includes: stent and disposable dilator

1 set per box, sterile

REQUIRES:

Disposable Sizer S1.1289u (see page 11)

S1.1390 25 mm

S1.1391 30 mm

S1.1392 35 mm

FCI
EXCLUSIVE



CRAWFORD


- For congenital nasolacrimal duct obstruction and canalicular pathologies
- Stainless steel probe with nontraumatic, rounded olive tip

Box of 3, sterile

REQUIRES:

Crawford Hook S1.1275 (see page 11)

S1.1270u



BIKA®


Designed by J.A. Bernard, M.D.

- For congenital nasolacrimal duct obstruction and canalicular pathologies
- Stainless steel probe with silicone tubing swaged into probe
- Nontraumatic rounded tip

Box of 3, sterile

S1.1000u

FCI
EXCLUSIVE



LACRIMAL STENTS & TUBES

Bicanalicular Intubation (cont.)

INFANT BIKA®


Designed by J.A. Bernard, M.D.

- For infant canalicular stenosis/DCR
- Stainless steel probes entirely embedded in silicone tubing
- Appropriate for DCR because of its high flexibility

Box of 3, sterile

S1.1100u

FCI
EXCLUSIVE



RITLENG®

Designed by P. Ritleng, M.D.

- For canalicular pathologies
- For nasal lacrimal duct obstruction
- Silicone, black PEEK (PolyEtherEtherKetone) thread-guide
- Silicone tubing (300 mm long, 0.64 mm diameter)
- Unique technique and design

Sterile

REQUIRES:

Ritleng Probe S1.1460 (see page 11)



See the surgical video and instructions at:
www.fci-ophthalmics.com/stents-tubes#ritleng

S1.1450u Box of 3

S1.1451u Box of 1

S1.1490u Box of 3 PVP

S1.1491u Box of 1 PVP

FCI
EXCLUSIVE



GOLDBERG CERCLAGE


Designed by R.A. Goldberg, M.D.

- For reconstruction, trauma and chronic stenosis of the canalicular system
- Silicone, blue polypropylene suture

Box of 3, sterile

S1.1900u

FCI
EXCLUSIVE



Monocanalicular Intubation

MASTERKA®

Designed by B. Fayet, M.D.

A unique combination of a punctum plug with a preloaded lacrimal duct stent.

- For congenital lacrimal duct obstructions
- Anchored at the punctum by a plug
- No intranasal recovery
- Silicone stent pre-mounted on stainless steel introducer

Sizing is important for proper stent selection


Box of 1, sterile

REQUIRES:

Disposable Punctum Dilator and Plug Inserter S1.3090u (see page 11)

Disposable Sizer S1.1289u (see page 11)

S1.1608	30 mm
S1.1609	35 mm
S1.1610	40 mm



MINI MONOKA®

Designed by B. Fayet, M.D. and J.A. Bernard, M.D.


- For canalicular lacerations involving external 2/3 of one canaliculus
- Stent is securely anchored at punctum by plug
- No knot or suture necessary
- Silicone

Box of 3, sterile

REQUIRES:

Disposable Punctum Dilator and Plug Inserter S1.3090u (see page 11)

S1.1500




WIDE/MEDIUM COLLARETTE MONOKA

Designed by B. Fayet, M.D. and J.A. Bernard, M.D.

- For congenital nasolacrimal duct obstruction and canalicular pathologies
- Silicone, stainless steel
- Nontraumatic rounded tip
- Silicone tubing swaged into malleable stainless steel probe
- Recommended: Disposable Punctum Dilator and Plug Inserter S1.3090u (see page 11)

Box of 3, sterile

S1.1630	4 mm	Wide
S1.1710	3 mm	Medium



LACRIMAL STENTS & TUBES

Monocanalicular Intubation (cont.)

SELF-THREADING MONOKA

Designed by B. Fayet, M.D., J.A. Bernard, M.D.
and P. Ritleng, M.D.

- For congenital nasolacrimal duct obstruction and canalicular pathologies
- Silicone, black PEEK (PolyEtherEtherKetone) thread-guide
- Securely anchored at punctum by plug
- No knot or suture necessary
- Recommended: Disposable Punctum Dilator and Plug Insertor S1.3090u (see page 11)


Sterile

REQUIRES:

Ritleng Probe S1.1460 (see page 11)

Ritleng Hook S1.1480 (see page 11)

S1.1810u	3 mm	Medium	Box of 3
S1.1811u	3 mm	Medium	Box of 1
S1.1800u	4 mm	Wide	Box of 3
S1.1801u	4 mm	Wide	Box of 1
S1.1820u	4 mm	Wide PVP	Box of 3
S1.1821u	4 mm	Wide PVP	Box of 1



MONO-CRAWFORD

- For congenital lacrimal duct obstruction
- For canalicular pathologies
- Single plug anchors stent at punctum
- Recommended: Disposable Punctum Dilator and Plug Insertor S1.3090u (page 11)

Sterile

REQUIRES:


Crawford Hook S1.1275 (see page 11)

3 MM

S1.1901u	Medium Collarette	Box of 1
S1.1903u	Medium Collarette	Box of 3
S1.1941u	PVP Medium Collarette	Box of 1
S1.1943u	PVP Medium Collarette	Box of 3

4 MM

S1.1911u	Wide Collarette	Box of 1
S1.1913u	Wide Collarette	Box of 3
S1.1951u	PVP Wide Collarette	Box of 1
S1.1953u	PVP Wide Collarette	Box of 3




Jones Tubes

STOPLOSS™

- Internal flexible silicone flange prevents tube extrusion
- Accurate tube sizing for each patient
- Simple to insert using dedicated disposable Introducer Set
- Large range of tube sizes
- Custom Removal Tool available

Sterile



FCI EXCLUSIVE

SLJT IS StopLoss™ Introducer Set

SLJT RT StopLoss™ Removal Tool

3.5 MM

SLJT 3509	9 mm
SLJT 3510	10 mm
SLJT 3511	11 mm
SLJT 3512	12 mm
SLJT 3513	13 mm
SLJT 3514	14 mm
SLJT 3515	15 mm
SLJT 3516	16 mm
SLJT 3517	17 mm
SLJT 3518	18 mm
SLJT 3519	19 mm
SLJT 3520	20 mm
SLJT 3521	21 mm
SLJT 3522	22 mm

4.0 MM

SLJT 4009	9 mm
SLJT 4010	10 mm
SLJT 4011	11 mm
SLJT 4012	12 mm
SLJT 4013	13 mm
SLJT 4014	14 mm
SLJT 4015	15 mm
SLJT 4016	16 mm
SLJT 4017	17 mm
SLJT 4018	18 mm
SLJT 4019	19 mm
SLJT 4020	20 mm
SLJT 4021	21 mm
SLJT 4022	22 mm

Lacrimal Stents & Tubes Accessories

RITLENG® PROBE

Box of 1, stainless steel, non-sterile

S1.1460 Ritleng Probe



ENDONASAL FORCEPS

Box of 1, stainless steel, non-sterile

S1.1470 Endonasal Forceps



FCI EXCLUSIVE

CRAWFORD HOOK

Box of 1, stainless steel, non-sterile

S1.1275 Crawford Hook



RITLENG® HOOK

Box of 1, stainless steel, non-sterile

S1.1480 Ritleng Hook



FCI EXCLUSIVE

DISPOSABLE SIZERS

Box of 5, stainless steel PolyOxy-Methylene (POM), sterile


S1.1289u Disposable Sizers



DISPOSABLE PUNCTUM DILATOR & PLUG INSERTER

Box of 3, sterile

S1.3090u Disposable Punctum Dilator & Plug Inserter



PIGTAIL PROBE

Box of 1, stainless steel, non-sterile

PT-1000 Pigtail Probe Standard Model

PT-2000 Pigtail Probe Pediatric Model



BOWMAN PROBE

Box of 1, stainless steel, non-sterile

BP-9-011 Bowman Probe Size 00 & 0

BP-9-012 Bowman Probe Size 1 & 2

BP-9-013 Bowman Probe Size 3 & 4



REUSABLE DOUBLE DILATOR

Box of 1, stainless steel, non-sterile

S2.3070 Fine & Medium Tapers



EYELID REPAIR

Ptosis Repair


PTOSE-UP

Designed by J.M. Ruban, M.D.

Ptose-Up is a suspension set consisting of a strip of biomaterial made of expanded Poly Tetra Fluoro Ethylene (ePTFE) fitted with a needle at each extremity. It is used to join the upper eyelid to the frontalis muscle in ptosis surgery with frontalis suspension.

- Strips made of a biocompatible, porous, inert, biointegratable, non-toxic, ready to use, non-allergenic material
- 58 mm long needles with ultra sharp tips to facilitate strip insertion
- Convenient to place and adjust
- Excellent eyelid contour
- Reduces surgery time
- Inert material that can be removed easily


Box of 1, sterile



S3.3021 2 mm x 150 mm

S3.3001 3 mm x 150 mm

FCI
EXCLUSIVE




PTOSIS SLING

Designed by J.L. Llouquet, M.D. and D. Guerin, M.D.

- For Ptosis (with frontalis suspension)
- Ultra-sharp stainless steel needles
- Needles: 58 mm long
- Smooth connection between silicone tubing and needles
- Silicone sleeve: outside diameter 1.65mm, inside diameter 0.76mm, 10mm in length

Box of 3, silicone sleeve in each blister, sterile

S3.1000u Ptosis Sling Set



Implants for Lagophthalmos

CONTOUR™ EYELID WEIGHTS

- Rapid recovery after operation and excellent tolerance
- Excellent functional and aesthetic results
- 1.0 mm thickness, 5.0 mm width
- Standard and special sizes
- Weight range from 0.6 grams to 2.8 grams

Box of 1, sterile

REQUIRED:

**Tantalum Eyelid Weight Sizing Set
TEW-7000 or TEW-7004 (see page 14)**

GOLD IMPLANTS

LL 3006	0.6 grams	7.6 mm Long
LL 3008	0.8 grams	10.2 mm Long
LL 3010	1.0 grams	12.3 mm Long
LL 3012	1.2 grams	14.5 mm Long
LL 3014	1.4 grams	16.6 mm Long
LL 3016	1.6 grams	18.8 mm Long
LL 3018	1.8 grams	21.0 mm Long
LL 3020	2.0 grams	23.1 mm Long
LL 3022	2.2 grams	25.2 mm Long
LL 3024	2.4 grams	27.3 mm Long
LL 3026	2.6 grams	29.4 mm Long
LL 3028	2.8 grams	31.5 mm Long



PLATINUM IMPLANTS

LL 5006	0.6 grams	7.6 mm Long
LL 5008	0.8 grams	10.2 mm Long
LL 5010	1.0 grams	12.3 mm Long
LL 5012	1.2 grams	14.5 mm Long
LL 5014	1.4 grams	16.6 mm Long
LL 5016	1.6 grams	18.8 mm Long
LL 5018	1.8 grams	21.0 mm Long
LL 5020	2.0 grams	23.1 mm Long
LL 5022	2.2 grams	25.2 mm Long
LL 5024	2.4 grams	27.3 mm Long
LL 5026	2.6 grams	29.4 mm Long
LL 5028	2.8 grams	31.5 mm Long



THINPROFILE™ EYELID WEIGHTS

- Implantable and external eyelid weights
- 40% thinner than Contour™ for more aesthetically pleasing treatment
- It provides a sleek, spherical curvature with tapered edges and rounded corners
- 0.65 mm thickness, weight from 0.6 grams to 1.8 grams
- Channel-free 1.0 mm holes for added freedom in suturing position
- Gold or platinum implants offer an excellent biocompatibility and functionality
- LL 4006-LL 4012 – 6.0 mm width
LL 4014-LL 4016 – 6.5 mm width
LL 6006-LL 6012 – 6.0 mm width
LL 6014-LL 6018 – 6.6 mm width

Box of 1, sterile

REQUIRED:

**Prior sizing. See Tantalum Eyelid Weight Sizing Sets
TEW-7000 or TEW-7004 (see page 14)**

GOLD IMPLANTS

LL 4006	0.6 grams	10.4 mm Long
LL 4008	0.8 grams	12.8 mm Long
LL 4010	1.0 grams	15.5 mm Long
LL 4012	1.2 grams	18.5 mm Long
LL 4014	1.4 grams	18.8 mm Long
LL 4016	1.6 grams	21.3 mm Long

PLATINUM IMPLANTS

LL 6006	0.6 grams	9.2 mm Long
LL 6008	0.8 grams	12.2 mm Long
LL 6010	1.0 grams	14.9 mm Long
LL 6012	1.2 grams	17.1 mm Long
LL 6014	1.4 grams	16.9 mm Long
LL 6016	1.6 grams	19.7 mm Long
LL 6018	1.8 grams	22.0 mm Long



EYELID REPAIR

External Weights for Lagophthalmos

BLINKEZE™ EYELID WEIGHTS

Includes Treatment Kit

- An immediate and easy way to treat temporary facial paralysis or Bell's palsy
- Can be quickly attached to the upper eyelid with the included double-sided and hypoallergenic adhesive strips (EyeClose® Adhesive Strips, LL 0505)
- After placement, the levator muscle seems stronger
- The proper weight allows to close the eyelid without difficulty and induces a slight ptosis (1.0 mm)
- Standard width 5.0 mm

Box of 1, non-sterile

REQUIRED:

Prior sizing and color selection.

See Tantalum Eyelid Weight Sizing Sets and Skin Tone Selector (see below)


**When placing an order, specify the corresponding color that best matches your patients' complexion (more information in FCI price list)*

STANDARD

TEW 6006	0.6 grams	17.1 mm Long
TEW 6008	0.8 grams	19.5 mm Long
TEW 6010	1.0 grams	11.9 mm Long
TEW 6012	1.2 grams	14.3 mm Long
TEW 6014	1.4 grams	16.7 mm Long
TEW 6016	1.6 grams	19.0 mm Long
TEW 6018	1.8 grams	21.4 mm Long

SPECIAL SIZES (CUSTOM ORDER*)

TEW 6020	2.0 grams	23.8 mm Long
TEW 6022	2.2 grams	26.2 mm Long
TEW 6024	2.4 grams	28.6 mm Long
TEW 6026	2.6 grams	30.9 mm Long
TEW 6028	2.8 grams	33.3 mm Long



TANTALUM EYELID WEIGHT SIZING SETS


- Sizing set to define the most accurate weight for the lid implant
- Recommended for all implants and external weights, including Blinkenze,™ Thinprofile™ and Coutour™

Set includes: one storage case, 5 or 7 weights, 1 skin tone selector card, 100 adhesive strips

TEW 7000 Standard Sizing Set, 7 Weights from 0.6 to 1.8 grams

TEW 7004 Special Sizing Set, 5 Weights from 2.0 to 2.8 grams

Skin Tone Selector



Eyelid Repair Accessories

OCULAR SHIELD

- Provided with suction cup for easy placement and removal
- Material: PMMA
- Single use

Box of 1, sterile

S6.2301



WRIGHT FASCIA NEEDLE

- Stainless steel
- Reusable

Box of 1, non-sterile

WF-1000u



Bioceramic Implants

The new standard in orbital implant material worldwide. These unique structures for vascularization and integration are bioinert, nontoxic, and nonallergenic.

MESH-WRAPPED IMPLANTS

Implant with a mesh to secure the attachment of extraocular muscles.


- For enucleation
- No need for donor sclera or autologous tissue
- Facilitated sphere introduction
- Mesh resorption in 45-60 days
- Material: aluminum oxide, vicryl mesh, silicone band
- Available in 4 diameters, from 16 mm to 22 mm

Box of 1, sterile, ready-to-use

MESH-WRAPPED

S6.5316u	16 mm
S6.5318u	18 mm
S6.5320u	20 mm
S6.5322u	22 mm

FCI
EXCLUSIVE



Mesh-Wrapped, Aluminum Oxide,
Vicryl Mesh, Silicone Band

SPHERICAL & EGG-SHAPED IMPLANTS

- For enucleation, evisceration or secondary implantation
- More biocompatible :
rapid and complete fibrovascularization
- Easy to handle and insert:
less sticky, no velcro effect
- Low tissue reaction and low exposure rate
- Porous, strong, non-brittle, stable, and lightweight
- 2 alternatives: sphere or egg-shaped
- 99% pure Aluminum Oxide (Al_2O_3)
- Porosity range: 200-550 μm

Box of 1, sterile

SPHERICAL

S6.5112u	12 mm
S6.5114u	14 mm
S6.5116u	16 mm
S6.5118u	18 mm
S6.5120u	20 mm
S6.5121u	21 mm
S6.5122u	22 mm

EGG-SHAPED

S6.5330u	18 mm
S6.5340u	20 mm

FCI
EXCLUSIVE



Spherical, Aluminum Oxide



Egg-Shaped, Aluminum Oxide

Silicone Implants

SPHERICAL IMPLANTS

- For enucleation or evisceration
- Medical grade silicone (translucent)
- 10 sizes, 12 mm to 22 mm

Box of 1, sterile

S6.1012u	12 mm
S6.1013u	13 mm
S6.1014u	14 mm
S6.1015u	15 mm
S6.1016u	16 mm
S6.1017u	17 mm
S6.1018u	18 mm
S6.1019u	19 mm
S6.1020u	20 mm
S6.1022u	22 mm



Spherical, Silicone

Self-Inflating Tissue Expanders

OSMED self-inflating tissue expanders are made of a specially developed hydrogel that uses the osmotic principle to gain volume.

These implants are indicated for the facilitation of normal growth of orbital tissues and bones to achieve proximate facial symmetry for the treatment of clinical anophthalmos, congenital microphthalmos, including post-enucleation socket syndrome.

After implantation, OSMED hydrogel implants start to absorb body fluid and grow consistently to a predefined form and volume.

The increased volume of the implant – between 6 to 12 fold, depending on the product type – leads to an increase of soft tissue.

- High biocompatibility
- Similar hydrogel also used in soft contact lenses
- Low risk of infection
- Controlled hydrogel swelling
- Does not use valves
- Does not require inflation or deflation

HEMISPHERE - SOCKET EXPANDER

For the treatment of clinical anophthalmia, congenial microphthalmia, and postenucleation implantation to achieve proximate facial symmetry.


- Expansion of the conjunctival sac with clinical anophthalmia
- Enlargement of the surface of the mucous membrane
- Creation of a fold abutment/support for artificial eye prosthesis

Box of 1, sterile

		Volume/Diameter Before Swelling	After Swelling*
OSM352-6040	0.4 ml	0.06 ml / 6 mm	0.4 ml / 11.2 mm
OSM352-6090	0.9 ml	0.13 ml / 8 mm	0.9 ml / 14.0 mm
OSM352-6090/P	0.9 ml/P	0.13 ml / 8 mm	0.9 ml / 14.0 mm
OSM352-6150	1.5 ml	0.20 ml / 9 mm	1.5 ml / 18.0 mm
OSM352-6200	2.0 ml	0.28 ml / 10 mm	2.0 ml / 20.0 mm

*in vitro in 0.9% NaCl-Sol.

FCI
EXCLUSIVE



SPHERE - ORBITA EXPANDER

For the treatment of clinical anophthalmia and to achieve proximate facial symmetry.


- Compensation of volume deficiency in the orbit
- Creation of an abutment/support for artificial eye prosthesis
- Enlargement of the palpebral fissure due to forward pressure from the artificial eye prosthesis

Box of 1, sterile

		Volume/Diameter Before Swelling	After Swelling*
OSM352-7100	1 ml	0.12 ml / 6 mm	1 ml / 12.4 mm
OSM352-7200	2 ml	0.30 ml / 8 mm	2 ml / 15.5 mm
OSM352-7300	3 ml	0.30 ml / 8 mm	3 ml / 18.0 mm
OSM352-7400	4 ml	0.43 ml / 9 mm	4 ml / 19.7 mm
OSM352-7500	5 ml	0.43 ml / 9 mm	5 ml / 21.8 mm

*in vitro in 0.9% NaCl-Sol.

FCI
EXCLUSIVE



Orbital Accessories

DISPOSABLE SIZERS

- For enucleation or evisceration
- To select accurate implant size
- PolyOxyMethylene (POM)
- Single use

Box of 4 in sizes 14, 16, 18 and 20 mm, sterile

S6.3060u



SPHERE INTRODUCER

- For inserting orbital implant into socket
- Stainless steel
- Reusable

Box of 1, sterile

S6.3050




SILICONE CONFORMER

- To maintain the conjunctival cul-de-sacs after enucleation or evisceration without hindering lid closure
- Holes allow drainage of secretions and administration of drops
- Silicone
- Translucent
- One size, can be trimmed
- Single use

Box of 1, sterile

S6.2001



Orbital Accessories (cont.)

PERFORATED CONFORMERS

- To maintain the conjunctival cul-de-sacs after enucleation or evisceration without hindering lid closure
- Holes allow drainage of secretions and administration of drops
- Polyamide 12
- Translucent for visual control of the conjunctiva
- Single use

Box of 1, sterile

S6.2230u	20 mm	Small
S6.2231u	23 mm	Medium
S6.2232u	26 mm	Large



SYMBLEPHARON RINGS

- To prevent retraction of the conjunctival cul-de-sacs without touching the cornea
- Edges are highly polished and atraumatic
- Polymethylmethacrylate (PMMA)
- Translucent for visual control of the conjunctiva
- Single use

Box of 1, sterile

S6.2420u	20 mm	Small
S6.2422u	22 mm	Medium
S6.2425u	25 mm	Large



Pre-Loaded CTRs in Disposable Injectors

EYEJET

- Choose between 2 insertion systems:
 - Right (for clockwise insertion)
 - Left (for counter clockwise insertion)
- Available for the following CTR types:
 - Morcher® 14, 14A, 14C
 - Henderson CTR
 - Malyugin/Cionni CTR



Visit [FCI-Ophthalmics.com](https://www.fci-ophthalmics.com) to watch the Eyejet informational video.

Capsular Tension Rings (CTR)

CTRs convert complex cataract cases from compromised zonular integrity into routine cases.

General Indications

- Missing or damaged zonules
- Lens subluxation
- Pseudoexfoliation
- Myopia
- Zonulolysis
- Soft IOL
- Marfan's Syndrome

General Advantages

- Circular expansion of capsular bags
- Stable conditions during surgery
- Improves IOL centration
- Reduced risk of capsular fibrosis
- Resists capsular shrinkage

MORCHER® CTR

- Designed to stabilize the capsular bag
- 3 different diameters to accommodate the variations in size of the capsule
- Unloaded, or pre-loaded on disposable injectors (Eyejets)
- Unloaded CTRs can be inserted manually or with the Geuder Injector G-32955 (see page 23)
- Opened and flexible filament of PMMA


Box of 1, sterile

UNLOADED CTR


MR-1400	Type 14
MR-1410	Type 14A
MR-1420	Type 14C

PRELOADED CTR ON INJECTOR (EYEJET)


MR-14 Right	Type 14
MR-14 Left	Type 14
MR-14A Right	Type 14A
MR-14A Left	Type 14A
MR-14C Right	Type 14C
MR-14C Left	Type 14C



Type 14
Axial length <24 mm
for normal eyes



Type 14A
Axial length >28 mm
for highly myopic eyes



Type 14C
Axial length 24-28 mm for normal
or myopic eyes

HENDERSON CTR

Designed by Bonnie Henderson, M.D.

Indentations allow access to cortical material that might be trapped by standard CTR.

- Designed for easier cortical removal while maintaining the desired stretch around the capsule
- Unloaded or pre-loaded on disposable injectors (Eyejets)
- Unloaded CTRs can be inserted manually or with the Geuder Injector G-32955 (see page 23)
- Opened and flexible filament of PMMA


Box of 1, sterile

UNLOADED CTR

MR-10C	Type 10C
---------------	----------

PRELOADED CTR ON INJECTOR (EYEJET)

MR-10C Right	Type 10C
---------------------	----------



Capsular Tension Rings (CTR) (cont.)

CIONNI CTR

Designed by Robert J. Cionni, M.D.

- Modified version of Morcher® CTR: specifically designed for scleral fixation with suture (eyelets)
- Cionni CTR used to stabilize and recenter the capsular bag
- Recommended: 9-0 Polypropylene Suture Non-Absorbable (see page 23)
- Opened and flexible filament of PMMA

Box of 1, sterile

Additional Indication

- More extreme zonular weakness or missing zonules larger than 4 clock hours

Additional Advantages

- May be sutured to the sclera
- Promotes circular expansion of the capsular bags


MR-1L Type 1L

MR-2C Type 2C*

MR-2L Type 2L

MR-1G Type 1G

*MR-2C may be inserted with Geuder Injector (see page 31), drawn up to 3/4 of the way only.



MALYUGIN/CIONNI CTR

Eliminates the risk of perforating the bag. The curved end allows to slide the CTR along the equator of the capsular bag.

- Eyelet at the curved end is sutured to the sclera for secure fixation and offers an alternative to capsule supporting hooks
- The only Cionni CTR preloaded on an Eyejet and in a clockwise injection mode
- Preloaded, 11 mm
- Recommended: 9-0 Polypropylene Suture Non-Absorbable (see page 23)
- Opened and flexible filament of PMMA

Box of 1, sterile

Additional Indication

- Severe or traumatic zonular dehiscence or dialysis with pupil abnormalities when Iris Hooks and other CTRs cannot achieve capsular stability

PRELOADED CTR ON INJECTOR (EYEJET)

MR-10G Right Type 10G



AHMED SEGMENT

Designed by Ike Ahmed, M.D.

- Arc of 120° available in two diameters and lengths
- Used to stabilize and recenter the capsular bag
- PMMA

Box of 1, sterile

Additional Indication

- Significant zonular weakness

Additional Advantages

- May be used intraoperatively (secured by iris retractor)
- May be used as permanent support (sutured to the sclera)

MR-6D	Type 6D	9.5 mm Diameter	9.61 mm Length
MR-6E	Type 6E	10.0 mm Diameter	10.14 mm Length



Cataract Accessories

GEUDER® CTR INJECTOR

- One-hand implantation by compressing the plunger
- Designed for clockwise insertion of Morcher® and Henderson CTRs
- Reusable

Box of 1, non-sterile

G-32955



GEUDER® PUPIL DILATOR INJECTOR

- Designed for clockwise one-hand implantation of the Morcher® Pupil Dilator (MR-5S)
- Reusable

Box of 1, non-sterile

G-32970



9-0 POLYPROPYLENE SUTURE NON-ABSORBABLE

- Stronger and particularly suitable for Cionni CTRs or for stabilizing an IOL
- Double-armed and packaged in box of 12
- Needles: straight or curved
- Suture: 5, 20 or 30 cm

16 MM STRAIGHT/16 MM 90° CURVED

P9009062 30 cm Suture

16 MM STRAIGHT/16 MM STRAIGHT

P90062 20 cm Suture

P90062-2 5 cm Suture

16 MM STRAIGHT/6.2 MM CURVED

P9006062 30 cm Suture

16 MM CURVED/16 MM CURVED

P99062 30 cm Suture



MORCHER® PUPIL DILATOR

- Semi-circular elastic ring for expansion of the pupil during phacoemulsification
- Optimal temporary implant for contracted pupils or floppy iris syndrome
- May be inserted manually or with the Geuder® Injector (see above)
- PMMA

Box of 1, sterile

MR-5S Type 5S

Overall Diameter: 7.50 mm

Pupil Size: 5.0–6.0 mm

Thickness: 0.15–0.60/0.90 mm



Important Note: The pupil dilator should be implanted only with viscoelastic for protection of the corneal endothelium.


Cataract Accessories (cont.)

RING CALIPERS

- Temporary template to create a more precise anterior capsulorhexis
- Polyethelmethacrilate (PEMA)

Box of 1, sterile

RC-5	5 mm Inner Diameter
RC-6	6 mm Inner Diameter



APX 200 PUPIL EXPANDER


APX 200 is a pupil expander that provides surgeons with an effective pupil dilation, comfortable surgical performance and a high safety profile, in eyes with small pupils or floppy iris (IFIS).

- Fast, easy, and safe positioning and removal
- Efficient iris expansion – quadrangular or trapezoidal shape according to the physician's preference and needs
- No intraocular manipulations required
- The physician does not operate over the device
- Does not interfere with surgical maneuvers
- Minimal learning curve
- One size, spring adjusts
- Nylon

Each set includes 2 forceps and 2 pupil expanders

APX 200-3 3 Sets per Box

APX 200-6 6 Sets per Box



IRIS RETRACTORS

Stabilizes the iris during cataract and retina surgery.

- Flexible blue polypropylene
- Adjustable silicone ring
- Reusable or single use

Sterile

DISPOSABLE

S9.5014 1 Set of 4

S9.5054 5 Sets of 4

REUSABLE

S9.5116 1 Set of 6, Packaged in Steam Autoclavable Container



MACKOOL CATARACT SUPPORT SYSTEM

Designed by Richard J. Mackool, M.D.

- Hooks provide capsular support during cataract surgery for patients with lax or absent zonules
- Comes with disposable stab knife
- Flexible and disposable
- Easy insertion and removal

Box of 5, sterile

MH-1005



Kitaro Kits

KITARO DRY/WET LAB

- Created with the learning surgeon in mind
- Simulates the stages of cataract surgery using both a wet and dry lab environment
- Ability to practice techniques repeatedly
- Removes the complications of long preparation times and costly live tissue samples
- Sold ready-to-use, however, optional goods such as a speculum and soft nucleus are available individually

Set includes a wide range of reusable and disposable goods such as the mask, forceps, anterior capsule film, and viscoelastic substance.

Contact FCI Ophthalmics for the full list of these items.

FCI
EXCLUSIVE

K-5000 KITARO Dry/Wet Lab



KITARO COMPLEX LABS

Two kits for advanced surgeons.

Lab1

Phakic IOL insertion or complex cases such as small pupil management, weak or ruptured zonules, and intraocular lens transcleral fixation.

Lab2

Conventional Extracapsular Extraction (ECCE) and posterior capsule rupture management.

Features and Benefits

- Ability to practice techniques and procedures with realistic synthetic materials
- Reusable and disposable goods are included in kits and sold individually as replacement parts only
- Sold ready-to-use, however, optional goods such as a drain box are available individually


Set includes a wide range of reusable and disposable goods such as the mask, forceps, corneas, and sclerocorneal sheets.

Contact FCI Ophthalmics for the full list of these items.

FCI
EXCLUSIVE

K-6000 KITARO Complex Lab1

K-7000 KITARO Complex Lab2




Autostable Trocar Cannula System

RETILOCK®

Unique Design for Optimum Stability in Complex Surgeries

- Exclusive and patented design to ensure stability of the cannula, with special scleral ring feature
- Ready-to-use for a one-step insertion
- Ultra sharp knife for smooth and high precision incision
- Cone-shaped cannula to make the insertion of instruments easier
- 3 mm and 4 mm scleral marker
- Supplied with a 23Ga and 25Ga infusion line compatible with BSS, gas and silicone oil



Trocar Cannula System

RETILOCK®

Valved

- 3 Trocar cannulas mounted on ultra sharp knives
- 1 infusion line
- Transparent silicone valve and cone-shaped cannula for easy insertion of instruments

Non-Valved

- 3 Trocar cannulas mounted on ultra sharp knives
- 1 infusion line
- Cone-shaped cannula for easy insertion of instruments
- 3 plugs: compatible with standard forceps such as Paufigue and Bonn forceps

Box of 5 kits, sterile, disposable

VALVED CANNULAS

S9.7100.23 23Ga

S9.7100.25 25Ga

NON-VALVED CANNULAS (Supplied with Plugs)

S9.7000.23 23Ga

S9.7000.25 25Ga

**FCI
EXCLUSIVE**



Valved



Non-Valved

*It is recommended to gently
insert the plug on the trocar.*

Probes

90° ADJUSTABLE LASER PROBES

Exclusive design

- Angle of curvature: pressure-controlled (from 0° to 90°)
- Safe probe removal: goes back in straight position when releasing the handle
- Compatible with most laser machines
- 50% stiffer than competition
- Short stroke for optimal control
- Low pressure needed to actuate the handle
- Ultra precise concentric laser spot
- Available in 23Ga, 25Ga

Box of 10, sterile

FCI
EXCLUSIVE

23GA

S9.1301.23	Iridex/Bausch & Lomb
S9.1302.23	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1303.23	Quantel
S9.1304.23	Nidek/Lumenis Novus Spectra
S9.1305.23	Dorc

25GA

S9.1301.25	Iridex/Bausch & Lomb
S9.1302.25	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1303.25	Quantel
S9.1304.25	Nidek/Lumenis Novus Spectra
S9.1305.25	Dorc



LASER PROBES

Ergonomic non-slip handpiece design

- Available in 20Ga, 23Ga, 25Ga, 27Ga
- Color coding for easy identification
- Ultra precise concentric laser spot
- Compatible with most laser machines

Box of 10, sterile



STRAIGHT TIP 20GA

S9.1001.20	Iridex/Bausch & Lomb
S9.1002.20	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1003.20	Quantel
S9.1004.20	Nidek/Lumenis Novus Spectra
S9.1005.20	Dorc

23GA

S9.1001.23	Iridex/Bausch & Lomb
S9.1002.23	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1003.23	Quantel
S9.1004.23	Nidek/Lumenis Novus Spectra
S9.1005.23	Dorc

25GA

S9.1001.25	Iridex/Bausch & Lomb
S9.1002.25	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1003.25	Quantel
S9.1004.25	Nidek/Lumenis Novus Spectra
S9.1005.25	Dorc

27GA

S9.1001.27	Iridex/Bausch & Lomb
S9.1002.27	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1003.27	Quantel
S9.1004.27	Nidek/Lumenis Novus Spectra
S9.1005.27	Dorc

CURVED TIP 20GA

S9.1101.20	Iridex/Bausch & Lomb
S9.1102.20	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1103.20	Quantel
S9.1104.20	Nidek/Lumenis Novus Spectra
S9.1105.20	Dorc

23GA

S9.1201.23	Iridex/Bausch & Lomb
S9.1202.23	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1203.23	Quantel
S9.1204.23	Nidek/Lumenis Novus Spectra
S9.1205.23	Dorc

25GA

S9.1201.25	Iridex/Bausch & Lomb
S9.1202.25	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1203.25	Quantel
S9.1204.25	Nidek/Lumenis Novus Spectra
S9.1205.25	Dorc

FLEXIBLE CURVED TIP 23GA

S9.1211.23	Iridex/Bausch & Lomb
S9.1212.23	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1213.23	Quantel
S9.1214.23	Nidek/Lumenis Novus Spectra
S9.1215.23	Dorc

25GA

S9.1211.25	Iridex/Bausch & Lomb
S9.1212.25	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1213.25	Quantel
S9.1214.25	Nidek/Lumenis Novus Spectra
S9.1215.25	Dorc

Probes (cont.)

LASER ILLUMINATION PROBES

Ergonomic non-slip handpiece design

- Available in 20Ga, 23Ga, 25Ga, 27Ga
- Color coding for easy identification
- Ultra precise concentric laser spot
- Compatible with most laser machines

Box of 10, sterile

19Ga Needle

- Unique design with universal 19Ga needle, to be inserted into FCI reusable light connector



Straight Illuminated Tip



Curved Illuminated Tip



19Ga Needle

STRAIGHT ILLUMINATED TIP 20GA

S9.1401.20	Iridex/Bausch & Lomb
S9.1402.20	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1403.20	Quantel
S9.1404.20	Nidek/Lumenis Novus Spectra
S9.1405.20	Dorc

23GA

S9.1401.23	Iridex/Bausch & Lomb
S9.1402.23	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1403.23	Quantel
S9.1404.23	Nidek/Lumenis Novus Spectra
S9.1405.23	Dorc

25GA

S9.1401.25	Iridex/Bausch & Lomb
S9.1402.25	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1403.25	Quantel
S9.1404.25	Nidek/Lumenis Novus Spectra
S9.1405.25	Dorc

27GA

S9.1401.27	Iridex/Bausch & Lomb
S9.1402.27	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1403.27	Quantel
S9.1404.27	Nidek/Lumenis Novus Spectra
S9.1405.27	Dorc

CURVED ILLUMINATED TIP 20GA

S9.1501.20	Iridex/Bausch & Lomb
S9.1502.20	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1503.20	Quantel
S9.1504.20	Nidek/Lumenis Novus Spectra
S9.1505.20	Dorc

23GA

S9.1601.23	Iridex/Bausch & Lomb
S9.1602.23	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1603.23	Quantel
S9.1604.23	Nidek/Lumenis Novus Spectra
S9.1605.23	Dorc

25GA

S9.1601.25	Iridex/Bausch & Lomb
S9.1602.25	Alcon/Coherent/Zeiss/Ellex/OS4
S9.1603.25	Quantel
S9.1604.25	Nidek/Lumenis Novus Spectra
S9.1605.25	Dorc

**All references include a universal 19Ga needle to fit with your light system*

Probes (cont.)

FIBER OPTIC PROBES

Ergonomic non-slip handpiece design


- Convenient 2.4 m long fibers
- Available in 20Ga, 23Ga, 25Ga, 27Ga

Focus, Widefield, Shielded Widefield: Box of 10, sterile

Chandelier: Box of 5, sterile

19Ga Needle

- Unique design with universal 19Ga needle, to be inserted into FCI reusable light connector



20GA

S9.2001.20	Focus Tip
S9.2101.20	Widefield Tip
S9.2201.20	Shielded Widefield Tip

23GA

S9.2001.23	Focus Tip
S9.2101.23	Widefield Tip
S9.2201.23	Shielded Widefield Tip

25GA

S9.2001.25	Focus Tip
S9.2101.25	Widefield Tip
S9.2201.25	Shielded Widefield Tip

25GA ALCON ACCURUS® OR CONSTELLATION®


S9.23.OE.25	Chandelier
--------------------	------------

25GA BAUSCH & LOMB STELLARIS PC®

S9.23.OD.25	Chandelier
--------------------	------------

27GA

S9.2001.27	Focus Tip
-------------------	-----------



LIGHT CONNECTORS

- Reusable

Box of 1, non-sterile

S9.24.0A	Alcon Acurus/Constellation®
S9.24.0B	Bausch & Lomb Millenium®
S9.24.0D	Bausch & Lomb Stellaris PC®
S9.24.0E	Dorc
S9.24.0F	Zeiss



Backflush Instruments

DISPOSABLE BACKFLUSH HANDPIECES

- All-in-one and ready-to-use
- Black silicone tips for better visibility
- Large reservoir for maximum reflux volume
- Plug for active or passive conversion
- Luer lock for active line attachment
- Available in 20Ga, 23Ga, 25Ga, 27Ga

Box of 5, sterile

20GA

S9.3012.20	Blunt Tip
S9.3010.20	Soft Tip
S9.3011.20	Brush Tip
S9.3013.20	Diamond Dusted Soft Tip

23GA

S9.3012.23	Blunt Tip
S9.3010.23	Soft Tip
S9.3011.23	Brush Tip
S9.3013.23	Diamond Dusted Soft Tip

25GA

S9.3012.25	Blunt Tip
S9.3010.25	Soft Tip
S9.3011.25	Brush Tip
S9.3013.25	Diamond Dusted Soft Tip

27GA

S9.3012.27	Blunt Tip
-------------------	-----------



Backflush with Cannula



Blunt Tip



Brush Tip



Soft Tip



Diamond Dusted Soft Tip

REUSABLE BACKFLUSH HANDLES

- Active or passive aspiration
- Durable titanium handle
- Designed to work with Disposable Reservoirs

Box of 1

REUSABLE HANDLES

S9.3110	Active Handle
S9.3120	Passive Handle

DISPOSABLE RESERVOIRS

S9.3210	Active Reservoir
S9.3220	Passive Reservoir



Active Handle



Passive Handle



Disposable Reservoirs

DISPOSABLE CANNULAS

- Several tip designs available

Box of 5, sterile

20GA

S9.4112.20	Blunt Tip
S9.4110.20	Soft Tip
S9.4111.20	Brush Tip
S9.4113.20	Dual Bore Tip

23GA

S9.4112.23	Blunt Tip
S9.4110.23	Soft Tip
S9.4111.23	Brush Tip
S9.4113.23	Dual Bore Tip
S9.4114.23	Subretinal Tip

25GA


S9.4112.25	Blunt Tip
S9.4110.25	Soft Tip
S9.4111.25	Brush Tip
S9.4113.25	Dual Bore Tip
S9.4014.25	Subretinal Tip




Blunt Tip



Soft Tip



Brush Tip



Dual Bore Tip



Subretinal Tip

Vitreoretinal Instruments

DISPOSABLE VITREORETINAL INSTRUMENTS

- Innovative anti-shredding tips design
- Maximum grasping power of 100 grams
- Low pressure needed to actuate the handle
- Short stroke for optimal control
- Exceptional precision: no forward or backward movement of the instrument tip
- Color coding for easy identification
- Available in 23Ga, 25Ga

Box of 5, sterile



FORCEPS 23GA

S9.6012.23 Eckardt (Maculorhexis)

S9.6013.23 End-Grasping

S9.6011.23 Serrated

S9.6017.23 Tano

25GA

S9.6012.25 Eckardt (Maculorhexis)

S9.6013.25 End-Grasping

S9.6011.25 Serrated

S9.6017.25 Tano

SCISSORS 23GA


S9.6014.23 Horizontal

S9.6015.23 Vertical


25GA

S9.6014.25 Horizontal


S9.6015.25 Vertical



Eckardt
(Maculorhexis)
Forceps



End-Grasping
Forceps




Serrated Forceps



Tano Forceps



Horizontal
Scissors



Vertical
Scissors

DISPOSABLE TANO POLISHER

- Color coding for easy identification
- Available in 20Ga, 23Ga, 25Ga
- Black diamonds and black silicone tip for better visibility
- Real diamond dust

Box of 5, sterile



20GA

S9.6016.20

23GA

S9.6016.23

25GA

S9.6016.25

Vitreoretinal Instruments (cont.)

REUSABLE VITREORETINAL INSTRUMENTS

Handles

- Made with lightweight and durable titanium
- Choice of one or two finger actuation
- Linear actuation for smooth closure

Tip Features

- Available in 17Ga, 20Ga, 23Ga, 25Ga
- Color coding for easy identification
- 360° rotating wheel for optimum control
- Supplied with Flushing Adapter A9.61.03

One Finger Control Handle

Two Finger Control Handle

Two Finger Control Handle with Tip

CONTROL HANDLES

A9.61.01 One Finger

A9.61.02 Two Finger

SCISSORS

20GA

A9.61.20.20 Straight, Rounded Blades

A9.61.21.20 Curved

A9.61.24.20 Angled 45° Horizontal

A9.61.26.20 Angled 55° Horizontal

A9.61.27.20 Angled 70° Vertical

A9.61.29.20 Klaus Lucke Retinotomy

23GA

A9.61.22.23 Curved

A9.61.27.23 Angled 70° Vertical

25GA

A9.61.22.25 Curved

A9.61.28.25 Angled 45° Vertical

MEMBRANE PICK

23GA

A9.61.43.23

25GA

A9.61.43.25

FORCEPS

17GA

A9.61.18.17 Foreign Body

20GA

A9.61.03.20 Diamonized, End-Gripping

A9.61.04.20 Tano Asymmetrical

A9.61.07.20 Serrated Jaws

A9.61.09.20 ILM

A9.61.08.20 Eckardt (Maculorhexis)

A9.61.16.20 Angled Horizontal, Diamonized

A9.61.05.20 Pick, Subretinal

23GA

A9.61.03.23 Diamonized, End-Gripping

A9.61.04.23 Tano Asymmetrical

A9.61.09.23 ILM

A9.61.08.23 Eckardt (Maculorhexis)

A9.61.05.23 Pick, Subretinal

25GA

A9.61.03.25 Diamonized, End-Gripping

A9.61.04.25 Tano Asymmetrical

A9.61.08.25 Eckardt (Maculorhexis)

A9.61.05.25 Pick, Subretinal



Membrane Pick



Foreign Body Forceps



Diamonized, End-Gripping Forceps



Tano Asymmetrical Forceps



Serrated Jaws Forceps



ILM Forceps



Eckardt (Maculorhexis) Forceps




Angled Horizontal, Diamonized Forceps



Pick Subretinal Forceps



Straight Scissors



Curved Scissors



Angled 45° Horizontal Scissors



Angled 55° Horizontal Scissors



Angled 70° Vertical Scissors



Klaus Lucke Retinotomy Scissors



Curved Scissors



Angled 45° Vertical Scissors

Scleral Buckling Components

SILICONE BANDS, STRIPS AND TIRES

FCI medical grade silicone implants for scleral buckling have been used in retinal detachment surgery throughout the world with great success.

- Soft and elastic
- Biocompatible and nontoxic
- Nonantigenic
- Easily tailored to exact size
- Wide range of shapes available

Sterile

BANDS (120 MM LENGTH)*

S5.1005u Type 40

S5.1015u Type 41

S5.1025u Type 42

S5.2005u Type 240

**All FCI bands packaged with fixation sleeves*

FIXATION SLEEVES*

S5.3005u Type 70 for Band Type 40

S5.3515u Type 71 for Band Type 41

S5.3025u Type 72 for Band Type 42

S5.3510u Type 270 for Band Type 240

**Supplied with bands and sold separately*

STRIPS (100 MM LENGTH)

S5.1105u Type 20

S5.1205u Type 31

S5.1255u Type 32

S5.2105u Type 219

S5.2205u Type 220

S5.2255u Type 225

TIRES

S5.2305u Type 277 90 mm

S5.2335u Type 279 110 mm

ASYMMETRICAL TIRES


S5.2295u Type 276 90 mm

S5.2345u Type 280 120 mm


BICONVEX TIRES

S5.2315u Type 287 90 mm


S5.2325u Type 289 120 mm




Type 40 Band




Type 41 Band




Type 42 Band




Type 240 Band




Type 70 Fixation Sleeve




Type 71 Fixation Sleeve




Type 72 Fixation Sleeve




Type 270 Fixation Sleeve




Type 20 Strip




Type 31 Strip




Type 32 Strip




Type 219 Strip




Type 220 Strip




Type 225 Strip




Type 276
Asymmetrical Tire




Type 277 Tire




Type 279 Tire



Type 280
Asymmetrical Tire



Type 287
Biconvex Tire



Type 289
Biconvex Tire

Scleral Buckling Components (cont.)

SILICONE SPONGES

FCI medical grade silicone implants for scleral buckling have been used in retinal detachment surgery throughout the world with great success.

- Soft and elastic
- Biocompatible and nontoxic
- Nonantigenic
- Easily tailored to exact size
- Wide range of shapes available

Sterile

SPONGE (80 MM LENGTH)

S5.6503u Type 503

S5.6504u Type 504

S5.6505u Type 505

S5.6535u Type 506

S5.6557u Type 507

S5.6585u Type 907

GROOVED SPONGE (80 MM LENGTH)

S5.6520u Type 505G

S5.6540u Type 506G

S5.6560u Type 507G

S5.6570u Type 509G


HALF SPONGE (80 MM LENGTH)

S5.6550u Type 510


S5.6575u Type 511

WATZKE SLEEVE SPREADING FORCEPS


S5.3100 Box of 1, Non-Sterile




Type 503 Sponge




Type 504 Sponge




Type 505 Sponge




Type 506 Sponge




Type 507 Sponge




Type 907 Sponge




Type 505G
Grooved Sponge




Type 506G
Grooved Sponge




Type 507G
Grooved Sponge



Type 509G
Grooved Sponge



Type 510
Half Sponge



Type 511
Half Sponge

Vitreoretinal Accessories

VITRECTOMY LENSES

Box of 5, sterile



Autostable Planoconcave



Planoconcave



Prismatic



Biconcave

AUTOSTABLE PMMA PLANOCONCAVE

S5.8050 30° x 0.94

PMMA

S5.8010u 30° x 0.94 Planoconcave

S5.8030u 32° x 0.94 35° Prismatic

S5.8040u 31.8° x 0.7 Biconcave

SILICONE

S5.7010u 30° x 1.08 Planoconcave

S5.7030u 30° x 1.08 30° Prismatic

S5.7040u 36° x 0.54 Biconcave

Vitreoretinal Accessories (cont.)

FCI 3-MIRROR LENS

Exceptional visualization every time

- Material: PMMA
- Angulation of 59°, 67°, 75°
- Central view 30°
- Image magnification 0.93x

Box of 5, sterile

\$5.8000



IRIS RETRACTORS

Stabilizes the iris during cataract and retina surgery.

- Flexible blue polypropylene
- Adjustable silicone ring
- Reusable or single use

Sterile

DISPOSABLE

- | | |
|----------------|-------------|
| S9.5014 | 1 Set of 4 |
| S9.5054 | 5 Sets of 4 |

REUSABLE

- | | |
|----------------|--|
| S9.5116 | 1 Set of 6, Packaged in Steam Autoclavable Container |
|----------------|--|



Intravitreal Injection Device

INVITRIA® INTRAVITREAL INJECTION ASSISTANT

Eliminates Unnecessary Instruments

No speculum, caliper or pressure plate needed.

Assures Accuracy

Fixed guide ensures proper angle, depth and distance of injection every time.

Relaxes Patients

InVitria® immobilizes the eye and hides the needle resulting in less patient anxiety.

Prevents Leakage

Rotating the device pre- and post-injection displaces and replaces the conjunctiva over the injection hole.

Compatible with:

- BD 30G 13 mm cannula
- B.Braun Sterican® 30G cannula


MW-200-I Box of 25, Sterile

Compatible with insulin-type cannulas and some safety needles such as the BD U-100 and the B.Braun 30G needles

MW-200-P Box of 25, Sterile

Compatible with BD 30G 13 mm cannula and B.Braun Sterican® 30G cannula

FCI EXCLUSIVE



DRY EYE

Permanent Punctum Plugs

Non-Sterile Plugs

Micro (wheel) S2.3524u	3
Mini (wheel) S2.3724u	3
Slim Mini (wheel) S2.3324u	3
Slim Petite (wheel) S2.3424u	3
Variety Pack (wheel) S2.3624u	3

PVP Perforated Plugs

Medium S1.3512u	4
Mini S1.3522u	4

Ready-Set Plugs

Large 1.0 mm S2.3522u	3
Medium 0.9 mm S2.3122u	3
Micro 0.6 mm S2.3652u	3
Mini 0.7 mm S2.3672u	3
Slim Mini 0.4 mm S2.3632u	3
Slim Petite 0.5 mm S2.3642u	3
Small 0.8 mm S2.3422u	3

SnugPlug™ S2.4002u

TearSaver Plugs

0.4 mm 814004	3
0.5 mm 814005	3
0.6 mm 814006	3
0.7 mm 814007	3
0.8 mm 814008	3

Temporary Punctum Plugs

Collagen Plugs

0.2 mm CP-0002	4
0.3 mm CP-0003	4
0.4 mm CP-0004	4
Variety Pack CP-0234	4

Prolong™ Absorbable Long Term Plugs

0.3 mm TP-3000	4
0.4 mm TP-4000	4
0.5 mm TP-5000	4
Variety Pack TP-3450	4

Tear Test

Zone Quick ZQ-1010	5
--------------------	---

Dry Eye Accessories

Punctal Gauges for Ready-Set Plugs S2.3060u	5
Reusable Double Dilator A8.4050	5

LACRIMAL STENTS & TUBES

Balloon Catheter

OphtaCath®

1 x 2 mm S1.4121	6
1 x 3 mm S1.4131	6
2 x 2 mm S1.4122	6
2 x 3 mm S1.4132	6

Bicanalicular Intubation

Bika® S1.1000u	7
Crawford S1.1270u	7
Goldberg Cerclage S1.1900	8
Infant Bika® S1.1100u	8
Nunchaku®	
Adult S1.1371	7
Children S1.1361	7

Bicanalicular Intubation (cont.)

Ritleng®

Box of 1 S1.1451u	8
Box of 3 S1.1450u	8
PVP - Box of 1 S1.1491u	8
PVP - Box of 3 S1.1490u	8

Self-Retaining Set II (SRS II)

25 mm S1.1390	7
30 mm S1.1391	7
35 mm S1.1392	7

Monocanalicular Intubation

Medium/Wide Collarette Monoka

Medium S1.1710	9
Wide S1.1630	9

Masterka®

30 mm S1.1608	9
35 mm S1.1609	9
40 mm S1.1610	9

Mini Monoka® S1.1500

Mono-Crawford/Medium Collarette

Box of 1 S1.1901u	10
Box of 3 S1.1903u	10
PVP - Box of 1 S1.1941u	10
PVP - Box of 3 S1.1943u	10

Mono-Crawford/Wide Collarette

Box of 1 S1.1911u	10
Box of 3 S1.1913u	10
PVP - Box of 1 S1.1951u	10
PVP - Box of 3 S1.1953u	10

Self-Threading Monoka/Medium

Box of 1 S1.1811u	10
Box of 3 S1.1810u	10

Self-Threading Monoka/Wide

Box of 1 S1.1801u	10
Box of 3 S1.1800u	10
PVP - Box of 1 S1.1821u	10
PVP - Box of 3 S1.1820u	10

Jones Tubes

StopLoss™ 9 mm - Flange

3.5 mm SLJT3509	10
4.0 mm SLJT4009	10

StopLoss™ 10 mm - Flange

3.5 mm SLJT3510	10
4.0 mm SLJT4010	10

StopLoss™ 11 mm - Flange

3.5 mm SLJT3511	10
4.0 mm SLJT4011	10

StopLoss™ 12 mm - Flange

3.5 mm SLJT3512	10
4.0 mm SLJT4012	10

StopLoss™ 13 mm - Flange

3.5 mm SLJT3513	10
4.0 mm SLJT4013	10

StopLoss™ 14 mm - Flange

3.5 mm SLJT3514	10
4.0 mm SLJT4014	10

StopLoss™ 15 mm - Flange

3.5 mm SLJT3515	10
4.0 mm SLJT4015	10

StopLoss™ 16 mm - Flange

3.5 mm SLJT3516	10
4.0 mm SLJT4016	10

StopLoss™ 17 mm - Flange

3.5 mm SLJT3517	10
4.0 mm SLJT4017	10

Jones Tubes (cont.)

StopLoss™ 18 mm - Flange	
3.5 mm SLJT3518	10
4.0 mm SLJT4018	10
StopLoss™ 19 mm - Flange	
3.5 mm SLJT3519	10
4.0 mm SLJT4019	10
StopLoss™ 20 mm - Flange	
3.5 mm SLJT3520	10
4.0 mm SLJT4020	10
StopLoss™ 21 mm - Flange	
3.5 mm SLJT3521	10
4.0 mm SLJT4021	10
StopLoss™ 22 mm - Flange	
3.5 mm SLJT3522	10
4.0 mm SLJT4022	10
StopLoss™ Introducer Set SLJT-IS	10
StopLoss™ Removal Tool SLJT-RT	10

Lacrimal Stents & Tubes Accessories

Bowman Probe	
Size 00 & 0 BP-9-011	11
Size 1 & 2 BP-9-012	11
Size 3 & 4 BP-9-013	11
Crawford Hook S1.1275	11
Disposable Punctum Dilator & Plug Inserter S1.3090u	11
Disposable Sizers (for SRS II & Masterka) S1.1289u	11
Endonasal Forceps S1.1470	11
Pigtail Probe	
Standard PT-1000	11
Pediatric PT-2000	11
Reusable Double Dilator S2.3070u	11
Ritleng Hook S1.1480	11
Ritleng Probe S1.1460	11

EYELID REPAIR

Ptosis Repair

Ptose-up	
2 mm x 150 mm S3.3021	12
3 mm x 150 mm S3.3001	12
Ptosis Sling S3.1000u	12

Implants for Lagophthalmos

Contour™ Eyelid Weights - Gold	
0.6 grams LL 3006	13
0.8 grams LL 3008	13
1.0 grams LL 3010	13
1.2 grams LL 3012	13
1.4 grams LL 3014	13
1.6 grams LL 3016	13
1.8 grams LL 3018	13
2.0 grams LL 3020	13
2.2 grams LL 3022	13
2.4 grams LL 3024	13
2.6 grams LL 3026	13
2.8 grams LL 3028	13
Contour™ Eyelid Weights - Platinum	
0.6 grams LL 5006	13
0.8 grams LL 5008	13
1.0 grams LL 5010	13
1.2 grams LL 5012	13
1.4 grams LL 5014	13
1.6 grams LL 5016	13

Implants for Lagophthalmos (cont.)

1.8 grams LL 5018	13
2.0 grams LL 5020	13
2.2 grams LL 5022	13
2.4 grams LL 5024	13
2.6 grams LL 5026	13
2.8 grams LL 5028	13
Thinprofile™ Eyelid Weights - Gold	
0.6 grams LL 4006	13
0.8 grams LL 4008	13
1.0 grams LL 4010	13
1.2 grams LL 4012	13
1.4 grams LL 4014	13
1.6 grams LL 4016	13
Thinprofile™ Eyelid Weights - Platinum	
0.6 grams LL 6006	13
0.8 grams LL 6008	13
1.0 grams LL 6010	13
1.2 grams LL 6012	13
1.4 grams LL 6014	13
1.6 grams LL 6016	13
1.8 grams LL 6018	13

External Weights for Lagophthalmos

Blinkeze™ Eyelid Weights	
0.6 grams TEW 6006	14
0.8 grams TEW 6008	14
1.0 grams TEW 6010	14
1.2 grams TEW 6012	14
1.4 grams TEW 6014	14
1.6 grams TEW 6016	14
1.8 grams TEW 6018	14
2.0 grams TEW 6020	14
2.2 grams TEW 6022	14
2.4 grams TEW 6024	14
2.6 grams TEW 6026	14
2.8 grams TEW 6028	14
EyeClose® Adhesive Strips LL 0505	14
Tantalum Eyelid Weight Sizing Set	
Special TEW 7004	14
Standard TEW 7000	14

Eyelid Repair Accessories

Ocular Shield S6.2301	15
Wright Fascia Needle WF-1000	15

ORBITAL

Bioceramic Implants

Mesh Wrapped Implants	
16 mm S6.5316u	16
18 mm S6.5318u	16
20 mm S6.5320u	16
22 mm S6.5322u	16
Egg-Shaped Implants	
18 mm S6.5330u	16
20 mm S6.5340u	16
Spherical Implants	
12 mm S6.5112u	16
14 mm S6.5114u	16
16 mm S6.5116u	16
18 mm S6.5118u	16
20 mm S6.5120u	16
21 mm S6.5121u	16
22 mm S6.5122u	16

Silicone Implants

Spherical Implants

12 mm S6.1012u	16
13 mm S6.1013u	16
14 mm S6.1014u	16
15 mm S6.1015u	16
16 mm S6.1016u	16
17 mm S6.1017u	16
18 mm S6.1018u	16
19 mm S6.1019u	16
20 mm S6.1020u	16
22 mm S6.1022u	16

Self-Inflating Tissue Expanders

Hemisphere

0.4 ml OSM352-6040	17
0.9 ml OSM352-6090	17
0.9 ml/P OSM352-6090/P	17
1.5 ml OSM352-6150	17
2.0 ml OSM352-6200	17

Sphere

1 ml OSM352-7100	17
2 ml OSM352-7200	17
3 ml OSM352-7300	17
4 ml OSM352-7400	17
5 ml OSM352-7500	17

Orbital Accessories

Disposable Sizers S6.3060u

Perforated Conformer

Large S6.2232u	19
Medium S6.2231u	19
Small S6.2230u	19

Silicone Conformer S6.2001

Sphere Introducer S6.3050

Symblepharon Ring

Large S6.2425u	19
Medium S6.2422u	19
Small S6.2420u	19

CATARACT

Capsular Tension Rings (CTR)

Ahmed Segment

Type 6D MR-6D	22
Type 6E MR-6E	22

Cionni

Type 1G MR-1G	22
Type 1L MR-1L	22
Type 2C MR-2C	22
Type 2L MR-2L	22

Morcher® Unloaded

Type 14 MR-1400	21
Type 14A MR-1410	21
Type 14C MR-1420	21

Morcher® Preloaded

Type 14 Right MR-14 Right	21
Type 14 Left MR-14 Left	21
Type 14A Right MR-14A Right	21
Type 14A Left MR-14A Left	21
Type 14C Right MR-14C Right	21
Type 14C Left MR-14C Left	21

Henderson Unloaded

Type 10C MR-10C	21
-----------------------	----

Henderson Preloaded

Type 10C Right MR-10C Right	21
-----------------------------------	----

Malyugin/Cionni Preloaded

Type 10G Right MR-10G Right	22
-----------------------------------	----

Cataract Accessories

9-0 Polypropylene Suture Non-Absorbable

Curved/Curved - 30 cm P99062	23
Straight/90° Curved - 30 cm P9009062	23
Straight/Curved - 30 cm P9006062	23
Straight/Straight - 5 cm P90062-2	23
Straight/Straight - 20 cm P90062	23

Morcher® Pupil Dilator MR-55

Ring Calipers

Type 5 RC-5	24
Type 6 RC-6	24

APX 200 Pupil Expander

3 Sets per Box APX 200-3	24
6 Sets per Box APX 200-6	24

Iris Retractors

Disposable - 1 Set of 4 S9.5014	24
Disposable - 5 Sets of 4 S9.5054	24
Reusable - 1 set of 6 S9.5116	24

Mackool Cataract Support System MH-1005

Geuder® CTR Injector G-32955

Geuder® Pupil Dilator Injector G-32970

Kitaro Kits

Kitaro Dry/Wet Lab K-5000	25
Kitaro Complex Lab1 K-6000	25
Kitaro Complex Lab2 K-7000	25

VITREORETINAL

Trocar Cannula System

Retilock®

23Ga - Non-Valved S9.7000.23	26
25Ga - Non-Valved S9.7000.25	26
23Ga - Valved S9.7100.23	26
25Ga - Valved S9.7100.25	26

Probes

90° Adjustable Laser Probes for Alcon/Coherent/Zeiss/Ellex/OS4

23Ga S9.1302.23	27
25Ga S9.1302.25	27

90° Adjustable Laser Probes for Dorc

23Ga S9.1305.23	27
25Ga S9.1305.25	27

90° Adjustable Laser Probes for Iridex/Bausch & Lomb

23Ga S9.1301.23	27
25Ga S9.1301.25	27

90° Adjustable Laser Probes for Nidek/Lumenis Novus Spectra

23Ga S9.1304.23	27
25Ga S9.1304.25	27

90° Adjustable Laser Probes for Quantel

23Ga S9.1303.23	27
25Ga S9.1303.25	27

Fiber Optic Probes - Chandelier

25Ga - Alcon Accurus® or Constellation® S9.23.OE.25	27
25Ga - Bausch & Lomb Stellaris PC® S9.23.OD.25	27

Fiber Optic Probes - Focus Tip

20Ga S9.2001.20	29
23Ga S9.2001.23	29
25Ga S9.2001.25	29
27Ga S9.2001.27	29

Fiber Optic Probe - Shielded Widefield Tip

20Ga S9.2201.20	29
23Ga S9.2201.23	29
25Ga S9.2201.25	29

Fiber Optic Probes - Widefield Tip

20Ga S9.2101.20	29
23Ga S9.2101.23	29
25Ga S9.2101.25	29

Probes (cont.)**Laser Probes for Alcon/Coherent/Zeiss/Ellex/OS4**

20Ga - Curved Tip S9.1102.20	27
20Ga - Straight Tip S9.1002.20	27
23Ga - Curved Tip S9.1202.23	27
23Ga - Flex Curved Tip S9.1212.23	27
23Ga - Straight Tip S9.1002.23	27
25Ga - Curved Tip S9.1202.25	27
25Ga - Flex Curved Tip S9.1212.25	27
25Ga - Straight Tip S9.1002.25	27
27Ga - Straight Tip S9.1002.27	27

Laser Probes for Dorc

20Ga - Curved Tip S9.1105.20	27
20Ga - Straight Tip S9.1005.20	27
23Ga - Curved Tip S9.1205.23	27
23Ga - Flex Curved Tip S9.1215.23	27
23Ga - Straight Tip S9.1005.23	27
25Ga - Curved Tip S9.1205.25	27
25Ga - Flex Curved Tip S9.1215.25	27
25Ga - Straight Tip S9.1005.25	27
27Ga - Straight Tip S9.1005.27	27

Laser Probes for Iridex/Bausch & Lomb

20Ga - Curved Tip S9.1101.20	27
20Ga - Straight Tip S9.1001.20	27
23Ga - Curved Tip S9.1201.23	27
23Ga - Flex Curved Tip S9.1211.23	27
23Ga - Straight Tip S9.1001.23	27
25Ga - Curved Tip S9.1201.25	27
25Ga - Flex Curved Tip S9.1211.25	27
25Ga - Straight Tip S9.1001.25	27
27Ga - Straight Tip S9.1001.27	27

Laser Probes for Nidek/Lumenis Novus Spectra

20Ga - Curved Tip S9.1104.20	27
20Ga - Straight Tip S9.1004.20	27
23Ga - Curved Tip S9.1204.23	27
23Ga - Flex Curved Tip S9.1214.23	27
23Ga - Straight Tip S9.1004.23	27
25Ga - Curved Tip S9.1204.25	27
25Ga - Flex Curved Tip S9.1214.25	27
25Ga - Straight Tip S9.1004.25	27
27Ga - Straight Tip S9.1004.27	27

Laser Probes for Quantel

20Ga - Curved Tip S9.1103.20	27
20Ga - Straight Tip S9.1003.20	27
23Ga - Curved Tip S9.1203.23	27
23Ga - Flex Curved Tip S9.1213.23	27
23Ga - Straight Tip S9.1003.23	27
25Ga - Curved Tip S9.1203.25	27
25Ga - Flex Curved Tip S9.1213.25	27
25Ga - Straight Tip S9.1003.25	27
27Ga - Straight Tip S9.1003.27	27

Laser Illumination Probes for Alcon/Coherent/Zeiss/Ellex/OS4

20Ga - Curved Tip S9.1502.20	28
20Ga - Straight Tip S9.1402.20	28
23Ga - Curved Tip S9.1602.23	28
23Ga - Straight Tip S9.1402.23	28
25Ga - Curved Tip S9.1602.25	28
25Ga - Straight Tip S9.1402.25	28
27Ga - Straight Tip S9.1402.27	28

Laser Illumination Probes for Dorc

20Ga - Curved Tip S9.1505.20	28
20Ga - Straight Tip S9.1405.20	28
23Ga - Curved Tip S9.1605.23	28
23Ga - Straight Tip S9.1405.23	28
25Ga - Curved Tip S9.1605.25	28
25Ga - Straight Tip S9.1405.25	28
27Ga - Straight Tip S9.1405.27	28

Probes (cont.)**Laser Illumination Probes for Iridex/Bausch & Lomb**

20Ga - Curved Tip S9.1501.20	28
20Ga - Straight Tip S9.1401.20	28
23Ga - Curved Tip S9.1601.23	28
23Ga - Straight Tip S9.1401.23	28
25Ga - Curved Tip S9.1601.25	28
25Ga - Straight Tip S9.1401.25	28
27Ga - Straight Tip S9.1401.27	28

Laser Illumination Probes for Nidek/Lumenis Novus Spectra

20Ga - Curved Tip S9.1504.20	28
20Ga - Straight Tip S9.1404.20	28
23Ga - Curved Tip S9.1604.23	28
23Ga - Straight Tip S9.1404.23	28
25Ga - Curved Tip S9.1604.25	28
25Ga - Straight Tip S9.1404.25	28
27Ga - Straight Tip S9.1404.27	28

Laser Illumination Probes for Quantel

20Ga - Curved Tip S9.1503.20	28
20Ga - Straight Tip S9.1403.20	28
23Ga - Curved Tip S9.1603.23	28
23Ga - Straight Tip S9.1403.23	28
25Ga - Curved Tip S9.1603.25	28
25Ga - Straight Tip S9.1403.25	28
27Ga - Straight Tip S9.1403.27	28

Light Connectors

Alcon Constellation® S9.24.0A	29
Bausch & Lomb Millenium® S9.24.0B	29
Bausch & Lomb Stellaris PC® S9.24.0D	29
Dorc S9.24.0E	29
Zeiss S9.24.0F	29

Backflush Instruments**Disposable Backflush Handpieces**

20Ga - Blunt Tip S9.3012.20	30
20Ga - Brush Tip S9.3011.20	30
20Ga - Diamond Dusted Soft Tip S9.3013.20	30
20Ga - Soft Tip S9.3010.20	30
23Ga - Blunt Tip S9.3012.23	30
23Ga - Brush Tip S9.3011.23	30
23Ga - Diamond Dusted Soft Tip S9.3013.23	30
23Ga - Soft Tip S9.3010.23	30
25Ga - Blunt Tip S9.3012.25	30
25Ga - Brush Tip S9.3011.25	30
25Ga - Diamond Dusted Soft Tip S9.3013.25	30
25Ga - Soft Tip S9.3010.25	30
27Ga - Blunt Tip S9.3012.27	30

Backflush Handles

Reusable Active Handle S9.3110	30
Reusable Passive Handle S9.3120	30

Disposable Reservoirs

Active Reservoir S9.3210	30
Passive Reservoir S9.3220	30

Disposable Cannulas

20Ga - Blunt Tip S9.4112.20	30
20Ga - Brush Tip S9.4111.20	30
20Ga - Dual Bore Tip S9.4113.20	30
20Ga - Soft Tip S9.4110.20	30
23Ga - Blunt Tip S9.4112.23	30
23Ga - Brush Tip S9.4111.23	30
23Ga - Dual Bore Tip S9.4113.23	30
23Ga - Soft Tip S9.4110.23	30
23Ga - Subretinal Tip S9.4114.23	30
25Ga - Blunt Tip S9.4112.25	30
25Ga - Brush Tip S9.4111.25	30
25Ga - Dual Bore Tip S9.4113.25	30
25Ga - Soft Tip S9.4110.25	30
25Ga - Subretinal Tip S9.4014.25	30

Vitreoretinal Instruments

Disposable Forceps	
23Ga - Eckardt (Maculorhexis) S9.6012.23	31
23Ga - End-Grasping S9.6013.23	31
23Ga - Serrated S9.6011.23	31
23Ga - Tano S9.6017.23	31
25Ga - Eckardt (Maculorhexis) S9.6012.25	31
25Ga - End-Grasping S9.6013.25	31
25Ga - Serrated S9.6011.25	31
25Ga - Tano S9.6017.25	31
Disposable Scissors	
23Ga - Horizontal S9.6014.23	31
23Ga - Vertical S9.6015.23	31
25Ga - Horizontal S9.6014.25	31
25Ga - Vertical S9.6015.25	31
Disposable Tano Polisher	
20Ga S9.6016.20	31
23Ga S9.6016.23	31
25Ga S9.6016.25	31
Reusable Control Handles	
One Finger A9.61.01	32
Two Finger A9.61.02	32
Reusable Membrane Pick	
23Ga A9.61.43.23	32
25Ga A9.61.43.25	32
Reusable Forceps	
17Ga - Foreign Body A9.61.18.17	32
20Ga - Angled Horizontal, Diamonized A9.61.16.20	32
20Ga - Diamonized, End-Gripping A9.61.03.20	32
20Ga - Eckardt (Maculorhexis) A9.61.08.20	32
20Ga - ILM A9.61.09.20	32
20Ga - Pick, Subretinal A9.61.05.20	32
20Ga - Serrated Jaws A9.61.07.20	32
20Ga - Tano Asymmetrical A9.61.04.20	32
23Ga - Diamonized, End-Gripping A9.61.03.23	32
23Ga - Eckardt (Maculorhexis) A9.61.08.23	32
23Ga - ILM A9.61.09.23	32
23Ga - Pick, Subretinal A9.61.05.23	32
23Ga - Tano Asymmetrical A9.61.04.23	32
25Ga - Diamonized, End-Gripping A9.61.03.25	32
25Ga - Eckardt (Maculorhexis) A9.61.08.25	32
25Ga - Pick, Subretinal A9.61.05.25	32
25Ga - Tano Asymmetrical A9.61.04.25	32
Reusable Scissors	
20Ga - Curved A9.61.21.20	32
20Ga - Angled 45° Horizontal A9.61.24.20	32
20Ga - Angled 55° Horizontal A9.61.26.20	32
20Ga - Klaus Lucke Retinotomy A9.61.29.20	32
20Ga - Straight, Rounded Blades A9.61.20.20	32
20Ga - Angled 70° Vertical A9.61.27.20	32
23Ga - Curved A9.61.22.23	32
23Ga - Angled 70° Vertical A9.61.27.23	32
25Ga - Curved A9.61.22.25	32
25Ga - Angled 45° Horizontal A9.61.28.25	32
Flushing Adapter A9.61.03	32

Scleral Buckling Components

Silicone Asymmetrical Tires	
Type 276 S5.2295u	33
Type 280 S5.2345u	33
Silicone Bands with Fixation Sleeve	
Type 40 S5.1005u	33
Type 41 S5.1015u	33
Type 42 S5.1025u	33
Type 240 S5.2005u	33

Scleral Buckling Components (cont.)

Silicone Biconvex Tires	
Type 287 S5.2315u	33
Type 289 S5.2325u	33
Silicone Fixation Sleeves	
Type 70 for Band Type 40 S5.3005u	33
Type 71 for Band Type 41 S5.3515u	33
Type 72 for Band Type 42 S5.3025u	33
Type 270 for Band Type 240 S5.3510u	33
Silicone Grooved Sponges	
Type 505G S5.6520u	34
Type 506G S5.6540u	34
Type 507G S5.6560u	34
Type 509G S5.6570u	34
Silicone Half Sponges	
Type 510 S5.6550u	34
Type 511 S5.6575u	34
Silicone Sponges	
Type 503 S5.6503u	34
Type 504 S5.6504u	34
Type 505 S5.6505u	34
Type 506 S5.6535u	34
Type 507 S5.6557u	34
Type 907 S5.6585u	34
Silicone Strips	
Type 20 S5.1105u	33
Type 31 S5.1205u	33
Type 32 S5.1255u	33
Type 219 S5.2105u	33
Type 220 S5.2205u	33
Type 225 S5.2255u	33
Silicone Tires	
Type 277 S5.2305u	33
Type 279 S5.2335u	33
Watzke Sleeve Spreading Forceps S5.3100	34
Vitreoretinal Accessories	
FCI 3-Mirror Lens S5.8000	35
Iris Retractors	
Disposable - 1 Set of 4 S9.5014	35
Disposable - 5 Sets of 4 S9.5054	35
Reusable - 1 set of 6 S9.5116	35
Vitrectomy Lenses	
Autostable PMMA Planoconcave S5.8050	34
PMMA Biconcave S5.8040u	34
PMMA Planoconcave S5.8010u	34
PMMA 35° Sign Prismatic S5.8030u	34
Silicone Biconcave S5.7040u	34
Silicone Planoconcave S5.7010u	34
Silicone 30° Sign Prismatic S5.7030u	34
Intravitreal Injection Device	
InVitria® Intravitreal Injection Assistant	
Compatible with Avastin® MW-200-I	35
Compatible with Lucentis® and Eylea® MW-200-P	35

All FCI products are completely LATEX-FREE.

FCI products are sold worldwide through a network of distributors and independent sales representatives. To place an order and/or to identify your distributor or representative, call toll-free 800.932.4202.

Payment Terms

30 days net

Shipping prepaid and added to invoice

Remit Payment to:

FCI Ophthalmics

30 Corporate Park Drive, Suite 310/320

Pembroke, MA 02359

accounts.receiveable@fci-ophthalmics.com

Returns

Defective merchandise may be returned for replacement only.

Merchandise shipped in error by FCI may be returned for credit or replacement within 30 days from date of invoice.

Note: Instruments and Blinkeze are non returnable.

- Prior authorization MUST be obtained.
- Returned merchandise must be properly packaged to prevent damage during shipping.
- Returned merchandise must be in original condition – not opened, marked, crushed, or otherwise damaged.

Restocking charge may apply.

Guarantee

FCI products are unconditionally guaranteed against manufacturing defects. Any product found to have a manufacturing defect will be replaced free of charge.

CONTACT FCI

Orders

orders@fci-ophthalmics.com

Telephone

800.932.4202

781.826.9060

Fax

781.826.9062

Information

info@fci-ophthalmics.com



FCI Ophthalmics

30 Corporate Park Drive, Suite 310/320

Pembroke, MA 02359

FCI-Ophthalmics.com

800.932.4202