

FOR IMMEDIATE RELEASE

NEWS

FCI OPHTHALMICS

P.O. BOX 465

Marshfield Hills, MA 02051

Tel: 800-932-4202 or 781-826-9060

Fax: 781-826-9062

info@fci-ophthalmics.com

www.fci-ophthalmics.com

Contact: Anne Bohsack

October, 2006


FCI Ophthalmics introduces The MONO-CRAWFORD For Monocanalicular Intubation

The Mono-Crawford is a unique combination of the successful Monoka and the familiar Crawford lacrimal intubation systems. One end of the stent features a silicone plug that anchors the system securely in place at the punctum. The other end is the olive tipped probe that facilitates retrieval from the nose using the Crawford hook (also available from FCI).

Thanks to the anchoring plug, the single tube which is cut in the nose requires no knots or sutures. Stent removal is easy – simply grasp the plug below it's collarette and pull out of the punctum.

The use of the Mono-Crawford is appropriate for congenital naso-lacrimal duct obstruction, laceration, or in cases where only one punctum is functional.

The Mono-Crawford stents come with either medium (3mm) or wide (4mm) collarettes. Stents are also available with a PVP coating which may enhance tear drainage while the stent is in place.



For more information contact:

Anne Bohsack @ 800-932-4202 or info@fci-ophthalmics.com

#