



# CATARACT TRAINING KIT

FCI-Ophthalmics.com  
30 Corporate Park Dr. #310/320  
Pembroke, MA 02359



## KITARO<sup>®</sup> EASY EYE

For IOL & CTR Insertion

K-4000	KITARO <sup>®</sup> Easy Eye (1-unit set)
K-4100	KITARO <sup>®</sup> Easy Eye (2-unit set)

### FEATURES AND BENEFITS

- Ability to practice techniques and procedures with realistic, synthetic materials
- For IOL and CTR practice
- The Easy Eye is composed of reusable and disposable goods
- Each set includes 1 movable stage, 1 sclera part, 1 sponge, 2 corneas and 2 posterior capsules
- Diameter: 50mm - Height: 48mm
- Additional corneas and capsules can be purchased separately

### PARTS YOU CAN USE WITH THE EASY EYE

- Corneas: K-2020 Cornea for WetLab and K-99 Cornea for iStent inject<sup>®</sup> (Glaucoma)
- Capsules: K-4010 Capsule for IOL and K-4010 Capsule for CTR insertion

To learn more about KITARO<sup>®</sup> Kits and accessories:

### Stay Connected



### Contact customer service at:

Phone: 800.932.4202

Fax: 781.826.9062

Email: [orders@fci-ophthalmics.com](mailto:orders@fci-ophthalmics.com)