



Morcher[®] EyeJet[®]

Available exclusively
from FCI

FCI-Ophthalmics.com
30 Corporate Park Dr. #310/320
Pembroke, MA 02359
800.932.4202



Preloaded Capsular Tension Rings

FEATURES AND BENEFITS

NEW UNIQUE INWARD BEND OF THE LEADING END, with all the features and benefits of the MORCHER[®] EYEJET[®] CTR Type 14.

The curved feature allows the CTR to slide along the equator of the capsular bag, thereby eliminating the risk of perforating the capsular fornix. Regular CTRs hit almost at an angle, but CTR Type 15, with one curved end helps:

- reduce possible puncturing (or bulging) of the capsular bag
- reduce risk of entrapment
- reduce risk of potentially tearing the posterior capsule
- reduce risk of inducing iatrogenic zonular damage during insertion

As it is being rotated into the bag, the EYEJET[®] CTR Type 15

- facilitates smooth insertion
- minimizes risk
- helps prevent complications

| Product # | Type | Equivalent Size | Compressed | Expanded | Axial Length | Eye | White to White |
|--------------|------|-----------------|------------|----------|--------------|-------------------|----------------|
| MR-15 Right | 15 | 14 | 10 mm | 12.3 | <24 mm | Normal | <11 mm |
| MR-15A Right | 15A | 14C | 11 mm | 13 | <24-28 mm | Normal/ myopic | 11-12.5 mm |
| MR-15B Right | 15B | 14A | 12 mm | 14.5 | >28 mm | Highly myopic | >12.5 mm |

Stay Connected



For more information contact customer service at:

Phone: 800.932.4202

Fax: 781.826.9062

Email: orders@fci-ophthalmics.com