



PRECIVIA[®]

THE GOLD STANDARD*
FOR ALL INTRAVITREAL INJECTIONS

FCI EXCLUSIVE



*Petri, A., Boysen, K., Cehofski, L.J. et al. Intravitreal Injections with Vascular Endothelial Growth Factor Inhibitors: A Practical Approach. Ophthalmol Ther 9, 191–203 (2020)

MORE THAN 1 MILLION INJECTIONS PERFORMED WORLDWIDE!

Precivia® simplifies intravitreal injections by making the procedure faster and more predictable. Its unique design presents a guide tube accurately oriented to ensure a fixed injection angle of 28°, a fixed distance at 3.5mm from the limbus, and a controlled depth.

LUCENTIS®¹ & EYLEA®²

- Lucentis® injection: prefilled syringe and vial format
- Eylea® injection: vial format
- Both compatible with BD 30Ga 13mm cannula and B. Braun Sterican™³ 30Ga cannula

		Single use
		Packaging
S5.200.P	Precivia®	Box of 25

AVASTIN®^{1*}

- Compatible with insulin type cannulas and safety needles such as BD U-100 and B. Braun 30Ga models

		Single use
		Packaging
S5.200.I	Precivia®	Box of 25

1. Lucentis® and Avastin® are registered trademarks of Genentech, Inc. 2. Eylea® is a registered trademark of Regeneron Pharmaceuticals, Inc. 3. Sterican™ is a trademark of B. Braun Melsungen AG. *Please verify the medical use of Avastin® with your local authorities.



FCI-Ophthalmics.com



FCI-Ophthalmics.com/contact