



NEW!

FCI-Ophthalmics.com
30 Corporate Park Dr. #310/320
Pembroke, MA 02359

SUPER GRIP FORCEPS

S9.6018.23 Super Grip - 23Ga

S9.6018.25 Super Grip - 25Ga

S9.6018.27 Super Grip - 27Ga

Innovative 3D Anti-Shredding Tip!

FCI Super Grip Forceps are specially designed to grasp the fibrous membranes that are difficult to catch with traditional instruments. A 360° actuation handle provides easy and precise maneuverability.

FEATURES AND BENEFITS

- Innovative 3D anti-shredding tip design with 100 grams grasping force
- Micro serrations on the broad platforms to grasp fibrous membranes
- Longer shaft for myopic eyes: 32mm length
- Low pressure needed to actuate 360° handle
- Removable handle for selective dexterity
- Short stroke for optimal control
- Available in 23, 25, and 27Ga
- Sterile, box of 5

To learn more about FCI's innovative vitreoretinal accessories:

Contact customer service at:

Phone: 800.932.4202

Fax: 781.826.9062

Email: orders@fci-ophthalmics.com



Stay Connected

