

# RetiLock®

## One-Step Self-Retaining Trocar System

**Non Valved & Valved**

For Easy Access to the  
Posterior Segment



800.932.4202

[www.FCI-Ophthalmics.com](http://www.FCI-Ophthalmics.com)

FCI OPHTHALMICS  
20 Winter Street • Pembroke, MA 02359

EMAIL  
[info@fci-ophthalmics.com](mailto:info@fci-ophthalmics.com)

TELEPHONE  
800.932.4202 • 781.826.9060

FAX  
781.826.9062



[www.FCI-Ophthalmics.com](http://www.FCI-Ophthalmics.com)

RLB0214



## SELF-RETAINING

- Exclusive design that ensures good stability of the cannula

## FEATURES & BENEFITS

- Ready-to-use one-step trocar system
- Mounted on an ultra sharp knife for smooth and high precision incision
- 3 and 4 mm scleral marker
- Supplied with a 23 and 25 Ga infusion line compatible with BSS, gas and silicone oil

## RETILOCK NON VALVED

### 23 & 25 GAUGE

- 3 trocar cannulas mounted on ultra sharp knives
- 1 infusion line
- Cone-shaped cannula for easy insertion of instruments
- 3 plugs: compatible with standard forceps such as Paufigue and Bonn Forceps. *It is recommended to gently insert the plug on the trocar.*



NON VALVED

## RETILOCK VALVED

### 23 & 25 GAUGE

- 3 trocar cannulas mounted on ultra sharp knives
- 1 infusion line
- Transparent silicone valve and cone-shaped cannula for easy insertion of instruments



VALVED

### NON VALVED 23 GAUGE

**S9-7000-23**  
Sterile  
Box of 5 kits



### NON VALVED 25 GAUGE

**S9-7000-25**  
Sterile  
Box of 5 kits



### VALVED 23 GAUGE

**S9-7100-23**  
Sterile  
Box of 5 kits



### VALVED 25 GAUGE

**S9-7100-25**  
Sterile  
Box of 5 kits

