

OCULID[®] Eyelid Weights

The Comprehensive Solution for Lagophthalmos

Effective temporary and permanent solutions for the management of functional defects and ocular exposure caused by blink, nocturnal, paralytic and iatrogenic lagophthalmos.



OCULID® IMPLANTABLE EYELID WEIGHTS

OCULID® Implantable Eyelid Weights conform to the surface curvature of the eye, offering superior cosmesis and aesthetic satisfaction among patients.

- Rapid post-operative recovery and excellent tolerance
- Excellent biocompatibility
- Restores dynamic blink and closure for functional results and lid cosmesis
- 1.0 mm thickness, 5.0 mm width
- Weight range from 0.6 g to 2.8 g
- 1.0 mm suture holes
- MRI Conditional**
- Sterile
- Box of 1



Required: Prior sizing.
See OCULID® Sizing Set.

Product Code	Material	Weight
S3.3006	OCULID Gold*	0.6 g
S3.3008	OCULID Gold*	0.8 g
S3.3010	OCULID Gold*	1.0 g
S3.3012	OCULID Gold*	1.2 g
S3.3014	OCULID Gold*	1.4 g
S3.3016	OCULID Gold*	1.6 g
S3.3018	OCULID Gold*	1.8 g
S3.3020*	OCULID Gold*	2.0 g
S3.3022*	OCULID Gold*	2.2 g
S3.3024*	OCULID Gold*	2.4 g
S3.3026*	OCULID Gold*	2.6 g
S3.3028*	OCULID Gold*	2.8 g

OCULID® SLIM IMPLANTABLE EYELID WEIGHTS

OCULID® Slim Implantable Eyelid Weights feature an anatomically friendly sleek surface with tapered edges and smoothly rounded corners. Available in gold and platinum.

- 40% thinner than OCULID® Standard Implantable Eyelid Weights
- Platinum density 11% higher than gold, providing smaller profile
- Weight range from 0.6 g to 1.8 g
- 1.0 mm suture holes
- MRI Conditional**
- Sterile
- Box of 1

Required: Prior sizing.
See OCULID® Sizing Set.

Product Code	Material	Weight
S3.4006	OCULID Gold Slim*	0.6 g
S3.4008	OCULID Gold Slim*	0.8 g
S3.4010	OCULID Gold Slim*	1.0 g
S3.4012	OCULID Gold Slim*	1.2 g
S3.4014	OCULID Gold Slim*	1.4 g
S3.4016	OCULID Gold Slim*	1.6 g
S3.6006	OCULID Platinum Slim*	0.6 g
S3.6008	OCULID Platinum Slim*	0.8 g
S3.6010	OCULID Platinum Slim*	1.0 g
S3.6012	OCULID Platinum Slim*	1.2 g
S3.6014	OCULID Platinum Slim*	1.4 g
S3.6016	OCULID Platinum Slim*	1.6 g
S3.6018	OCULID Platinum Slim*	1.8 g



*Special sizes available through custom order

**Check IFU for MRI Parameters

OCULID® SIZING SET

OCULID® Sizing Set assists in selecting the appropriate implantable or external eyelid weight. The eyelid should close without difficulty and induce a slight ptosis of 1.0 mm.

Set includes:

- A storage case
- 7 sizing weights from 0.6 g to 1.8 g (.2 g increments)
- 100 double-sided, hypoallergenic adhesive strips
- Skin tone selector card to best match patient's complexion
- For external use only
- Reusable, non-sterile

Product Code	Packaging
S3.7000	Box of 1



HYPALLERGENIC ADHESIVE STRIPS

To attach OCULID® External Eyelid Weights to the upper eyelid.

Product Code	Packaging
S3.0500	Box of 100



OCULID® EXTERNAL EYELID WEIGHTS

- Effective non-surgical treatment
- Simple use and application
- Allows eyelid closure without difficulty
- Hidden behind superior palpebral fold
- Helps protect surface of the eye
- Relieves dry eye
- Strengthens levator muscles with continued use
- Induces slight ptosis (1.0 mm)
- Non-sterile
- Box of 1

Required: Prior sizing. See OCULID® Sizing Set.

Product Code	Color	Weight
S3.6006.0	Light Pink	0.6 g
S3.6008.0	Light Pink	0.8 g
S3.6010.0	Light Pink	1.0 g
S3.6012.0	Light Pink	1.2 g
S3.6014.0	Light Pink	1.4 g
S3.6016.0	Light Pink	1.6 g
S3.6018.0	Light Pink	1.8 g
S3.6006.1	Pink-Neutral	0.6 g
S3.6008.1	Pink-Neutral	0.8 g
S3.6010.1	Pink-Neutral	1.0 g
S3.6012.1	Pink-Neutral	1.2 g
S3.6014.1	Pink-Neutral	1.4 g
S3.6016.1	Pink-Neutral	1.6 g
S3.6018.1	Pink-Neutral	1.8 g
S3.6006.2	Tan	0.6 g
S3.6008.2	Tan	0.8 g
S3.6010.2	Tan	1.0 g
S3.6012.2	Tan	1.2 g
S3.6014.2	Tan	1.4 g
S3.6016.2	Tan	1.6 g
S3.6018.2	Tan	1.8 g

Product Code	Color	Weight
S3.6006.3	Medium Brown	0.6 g
S3.6008.3	Medium Brown	0.8 g
S3.6010.3	Medium Brown	1.0 g
S3.6012.3	Medium Brown	1.2 g
S3.6014.3	Medium Brown	1.4 g
S3.6016.3	Medium Brown	1.6 g
S3.6018.3	Medium Brown	1.8 g
S3.6006.4	Dark Brown	0.6 g
S3.6008.4	Dark Brown	0.8 g
S3.6010.4	Dark Brown	1.0 g
S3.6012.4	Dark Brown	1.2 g
S3.6014.4	Dark Brown	1.4 g
S3.6016.4	Dark Brown	1.6 g
S3.6018.4	Dark Brown	1.8 g

Stay Connected!



Social Media



Learn More!



FCI®

a ZEISS company

FCI OPHTHALMICS

30 Corporate Park Drive, Suite 310/320
Pembroke, MA 02359

EMAIL

info@fci-ophthalmics.com

TELEPHONE

800.932.4202 • 781.826.9060

FAX

781.826.9062

www.FCI-Ophthalmics.com