

OCULID[®] Eyelid Weights

The Comprehensive Solution for Lagophthalmos

NEW!

Effective temporary and permanent solutions for the management of functional defects and ocular exposure caused by blink, nocturnal, paralytic and iatrogenic lagophthalmos.



www.FCI-Ophthalmics.com



OCULID GOLD® IMPLANTS

OCULID® implants conform to the surface curvature of the eye, offering superior cosmesis and aesthetic satisfaction among patients. Additionally, they feature a sleek surface, anatomically friendly curvature with tapered edges, and smoothly rounded corners.

- Rapid post-operative recovery and excellent tolerance
- Excellent biocompatibility
- Restores dynamic blink and closure for functional results and lid cosmesis
- 1.0 mm thickness, 5.0 mm width
- Weight range from 0.6 g to 2.8 g in 0.2 g increments*
- 1.0 mm suture holes
- MRI Conditional**
- Sterile
- Box of 1



Product Code	Material	Weight
S3.3006	Gold	0.6 g
S3.3008	Gold	0.8 g
S3.3010	Gold	1.0 g
S3.3012	Gold	1.2 g
S3.3014	Gold	1.4 g
S3.3016	Gold	1.6 g
S3.3018	Gold	1.8 g
S3.3020*	Gold	2.0 g
S3.3022*	Gold	2.2 g
S3.3024*	Gold	2.4 g
S3.3026*	Gold	2.6 g
S3.3028*	Gold	2.8 g

Required: Prior sizing. See OCULID® Sizing Set.

*Sizes 2.0 g to 2.8 g are available through custom order

**Check IFU for MRI Parameters

OCULID GOLD SLIM® | OCULID PLATINUM SLIM® IMPLANTS

OCULID® implants conform to the surface curvature of the eye, offering superior cosmesis and aesthetic satisfaction among patients. Additionally, they feature a sleek surface, anatomically friendly curvature with tapered edges, and smoothly rounded corners.

- Rapid post-operative recovery and excellent tolerance
- Excellent biocompatibility
- Restores dynamic blink and closure for functional results and lid cosmesis
- 0.6 mm thickness, 5.0 mm width
- Weight range from 0.6 g to 1.8 g
- Available in gold and platinum
- 40% thinner than OCULID Gold® implants
- Platinum density 11% higher than gold
- 1.0 mm suture holes
- MRI Conditional**
- Sterile
- Box of 1



Product Code	Material	Weight
S3.4006	Gold	0.6 g
S3.4008	Gold	0.8 g
S3.4010	Gold	1.0 g
S3.4012	Gold	1.2 g
S3.4014	Gold	1.4 g
S3.4016	Gold	1.6 g
S3.6006	Platinum	0.6 g
S3.6008	Platinum	0.8 g
S3.6010	Platinum	1.0 g
S3.6012	Platinum	1.2 g
S3.6014	Platinum	1.4 g
S3.6016	Platinum	1.6 g
S3.6018	Platinum	1.8 g

Required: Prior sizing. See OCULID® Sizing Set.

**Check IFU for MRI Parameters

OCULID® SIZING SET

OCULID® Sizing Set assists in selecting the appropriate implantable or external eyelid weight. The eyelid should close without difficulty and induce a slight ptosis of 1.0 mm.

Set includes:

- A storage case
- 7 sizing weights from 0.6 g to 1.8 g in 0.2 g increments
- 100 double-sided, hypoallergenic adhesive strips
- Skin tone selector card to best match patient's complexion
- For external use only
- Reusable, non-sterile

Product Code	Packaging
S3.7000	Box of 1



HYPOALLERGENIC ADHESIVE STRIPS

To attach OCULID® External Eyelid Weights to the upper eyelid.



Product Code	Packaging
S3.0500	Box of 100

OCULID® EXTERNAL EYELID WEIGHTS

- Effective non-surgical treatment
- Simple use and application
- Allows eyelid closure without difficulty
- Hidden behind superior palpebral fold
- Helps protect surface of the eye
- Relieves dry eye
- Strengthens levator muscles with continued use
- Induces slight ptosis (1.0 mm)
- Non-sterile
- Box of 1

Required: Prior sizing. See OCULID® Sizing Set.


Product Code	Color	Weight
S3.6006-0	Light Pink	0.6 g
S3.6008-0	Light Pink	0.8 g
S3.6010-0	Light Pink	1.0 g
S3.6012-0	Light Pink	1.2 g
S3.6014-0	Light Pink	1.4 g
S3.6016-0	Light Pink	1.6 g
S3.6018-0	Light Pink	1.8 g
S3.6006-1	Pink-Neutral	0.6 g
S3.6008-1	Pink-Neutral	0.8 g
S3.6010-1	Pink-Neutral	1.0 g
S3.6012-1	Pink-Neutral	1.2 g
S3.6014-1	Pink-Neutral	1.4 g
S3.6016-1	Pink-Neutral	1.6 g
S3.6018-1	Pink-Neutral	1.8 g
S3.6006-2	Tan	0.6 g
S3.6008-2	Tan	0.8 g
S3.6010-2	Tan	1.0 g
S3.6012-2	Tan	1.2 g
S3.6014-2	Tan	1.4 g
S3.6016-2	Tan	1.6 g
S3.6018-2	Tan	1.8 g

Product Code	Color	Weight
S3.6006-3	Medium Brown	0.6 g
S3.6008-3	Medium Brown	0.8 g
S3.6010-3	Medium Brown	1.0 g
S3.6012-3	Medium Brown	1.2 g
S3.6014-3	Medium Brown	1.4 g
S3.6016-3	Medium Brown	1.6 g
S3.6018-3	Medium Brown	1.8 g
S3.6006-4	Dark Brown	0.6 g
S3.6008-4	Dark Brown	0.8 g
S3.6010-4	Dark Brown	1.0 g
S3.6012-4	Dark Brown	1.2 g
S3.6014-4	Dark Brown	1.4 g
S3.6016-4	Dark Brown	1.6 g
S3.6018-4	Dark Brown	1.8 g

Stay Connected!



Social Media



Learn More!



FCI OPHTHALMICS

30 Corporate Park Drive, Suite 310/320
Pembroke, MA 02359

EMAIL

info@fci-ophthalmics.com

TELEPHONE

800.932.4202 • 781.826.9060

FAX

781.826.9062

www.FCI-Ophthalmics.com