


Punctal Plugs

The Comprehensive Solution for Dry Eye

INDICATIONS

Temporary and Permanent Punctal Plugs
for Dry Eye Management



FCI®

a ZEISS company



www.FCI-Ophthalmics.com

TEMPORARY PUNCTAL PLUGS

FCI PROLONG® DX

- Indicated for use in the diagnosis of dry eye
- To assess the effectiveness of permanent punctal occlusion
- Dissolves within **5 days**
- 2 mm length
- Sterile
- Made of collagen

6 plugs per packet, 10 packets per box (60 total)

Variety Pack: 6 plugs per packet–3 packets of 0.2 mm, 3 packets 0.3 mm, 4 packets of 0.4 mm (60 total)



S2.5202	0.2 mm	10 Packets of 6
S2.5203	0.3 mm	10 Packets of 6
S2.5204	0.4 mm	10 Packets of 6
S2.5224	Variety Pack	10 Packets of 6



FCI PROLONG®

- Indicated for use in the treatment of dry eye syndrome and dry eye conditions due to various ocular surface diseases
- For post-surgical dry eye symptoms
- Absorbs within **90 days**
- 2 mm length
- Sterile
- Made of copolymer

2 plugs per packet, 10 packets per box (20 total)

Variety Pack: 2 plugs per packet–3 packets 0.3 mm, 5 packets 0.4 mm, 2 packets 0.5 mm (20 total)



S2.5002	0.2 mm	10 Packets of 2
S2.5003	0.3 mm	10 Packets of 2
S2.5004	0.4 mm	10 Packets of 2
S2.5005	0.5 mm	10 Packets of 2
S2.5035	Variety Pack	10 Packets of 2



FCI PROLONG® +

- Indicated for use in the treatment of dry eye syndrome as a result of seasonal allergies, ocular surgery, or contact lens wear
- Absorbs within **180 days**
- 2 mm length
- Sterile
- Made of E-Caprolactone-L-Lactide copolymer (PCL)

2 plugs per packet, 10 packets per box (20 total)

Variety Pack: 2 plugs per packet–3 packets 0.3 mm, 4 packets 0.4 mm, 3 packets 0.5 mm (20 total)



S2.5103	0.3 mm	10 Packets of 2
S2.5104	0.4 mm	10 Packets of 2
S2.5105	0.5 mm	10 Packets of 2
S2.5135	Variety Pack	10 Packets of 2



PERMANENT PUNCTAL PLUGS

Designed for ease of use and indicated for punctal occlusion, FCI permanent punctal plugs are preloaded on a single-use inserter/dilator to facilitate insertion.

SNUGPLUG®

- One-size-fits-most
- No sizing needed
- High retention rate
- No dilation needed
- Ultra-thin collarette
- Medical grade silicone
- Individually packaged, sterile

S2.4002u One-size-fits-most Box of 2



Pre-stretched SNUGPLUG®

PRELOADED PUNCTAL PLUGS

DESIGNED BY J.A. BERNARD, M.D.

- Ready-set plugs
- Original design
- Ergonomic slanted collarette (except Slim Plugs)
- Good tolerance and optimal patient comfort
- Medical grade silicone
- Individually packaged, sterile

S2.3632	Slim Mini	0.4 mm	Box of 2
S2.3642	Slim Petite	0.5 mm	Box of 2
S2.3652	Extra Small	0.6 mm	Box of 2
S2.3672	Mini	0.7 mm	Box of 2
S2.3422	Small	0.8 mm	Box of 2
S2.3122	Regular	0.9 mm	Box of 2
S2.3522	Large	1.0 mm	Box of 2



The stamped likeness of the plug on the inserter indicates the direction of the collarette slant.



TEARSAVER PLUGS®

- Tapered shaft design
- Secure anchoring of the plug in the punctum
- More comfort and less risk of rub-out
- Permanent yet removable
- Medical grade silicone
- Individually packaged, sterile

S2.3842	0.4 mm	Box of 2
S2.3852	0.5 mm	Box of 2
S2.3862	0.6 mm	Box of 2
S2.3872	0.7 mm	Box of 2
S2.3882	0.8 mm	Box of 2




PERMANENT PUNCTAL PLUGS

NON-STERILE PLUGS

- Indicated for punctal occlusion
- Pre-mounted on a wheel for quick and easy loading onto a plug inserter
- Includes 2 Single Use Punctum Dilator and Plug Inserters (\$1.3090)
- Medical grade silicone
- Non-sterile

S2.3324u	Slim Mini	0.4 mm	Box of 24
S2.3424u	Slim Petite	0.5 mm	Box of 24
S2.3524u	Extra Small	0.6 mm	Box of 24
S2.3724u	Mini	0.7 mm	Box of 24
S2.3624u	6 plugs of each size	0.4, 0.5, 0.6, 0.7 mm	Box of 24



Variety Pack

DRY EYE ACCESSORIES

REUSABLE DOUBLE DILATOR

- Stainless steel
- Non-sterile

A8.4050	Fine & Medium	Box of 1
----------------	---------------	----------



STAY CONNECTED!



Social Media



Learn More!



FCI®
a ZEISS company

FCI OPHTHALMICS
30 Corporate Park Drive, Suite 310/320
Pembroke, MA 02359

EMAIL
info@fci-ophthalmics.com

TELEPHONE
800.932.4202 • 781.826.9060

FAX
781.826.9062

www.FCI-Ophthalmics.com